

CRL for Cereals and Feeding stuff
National Food Institute
Technical University of Denmark



Validation Report 6

Determination of pesticide residues in cereals by GC-MS/MS and LC-MS/MS

(QuEChERS method)

**Karen Hjorth,
Mette Erecius Poulsen &
Hanne Bjerre Christensen
10 March 2010**

CONTENT:

<i>1. Introduction</i>	2
<i>2. Principle of analysis</i>	2
<i>3. Validation design</i>	3
<i>4. Chromatograms and calibration curves</i>	4
<i>5. Validation parameters</i>	6
<i>6. Criteria for the acceptance of validation results</i>	6
<i>7. Results and discussion</i>	7
<i>8. Conclusions</i>	10
<i>9. References</i>	10
<i>Appendix 1. MRM transitions for the pesticides sought validated.</i>	11
<i>Appendix 2. Repeatability, recovery and limit of quantification for compounds validated on LC-MS-MS.</i>	15
<i>Appendix 3. Repeatability, recovery and limit of quantification for compounds validated on GC-MS-MS.</i>	18

1. Introduction

This report describes the validation of the QuEChERS method combined with GC-MS/MS and LC-MS/MS. The method was sought validated for 55 pesticides, isomers and degradation products in wheat, rye and oat. The method has previously been validated for about 98 pesticides.

The QuEChERS method has an extraction and clean-up step, which has been developed to be Quick, Easy, Cheap, Efficient, Rugged and Safe. The method is most commonly used on fruit and vegetables¹.

2. Principle of analysis

Sample preparation:

Cold water/ice water and acetonitril are added to the milled sample.

Extraction:

The sample is shaken and a salt and buffer mixture is added and the sample is shaken again.

Clean-up:

After centrifugation the supernatant is transferred to a tube with PSA and MgSO₄. After shaking and an additional centrifugation step the final extract is obtained.

For the LC-MS/MS analysis the extraction is followed by adding internal standard and the extract is filtered into HPLC vials.

Quantification and qualification:

GC-MS/MS:

The final extract is analysed GC/MS/MS (electron energy 70 eV, source temp. 180°C, transfer line GC interface 250°C). The injection volume is 4 µl.

LC-MS/MS:

The pesticide residues are separated on a reversed-phase column and detected by tandem mass spectrometry (MS/MS) by electrospray (ESI).

Selectivity and specificity:

GC-MS/MS:

GC-MS/MS is a highly selective method, and thereby highly specific. Two MRM transitions were used (two parent and two daughter ion) one for quantification and another transition for qualification. Parent and daughter ions are presented in appendix 1.

LC-MS/MS:

The validation includes pesticides determined with both positive and negative ESI. ¹³C₆-carbaryl was used as internal standard for quantification. All pesticides were detected in the multiple reaction monitoring mode (MRM). For each pesticide precursor ion and 2 product ions (where possible) were determined. One product ion for quantification and one for qualification. The MRM transitions for the pesticides and degradation products sought validated are given in appendix 1.

3. Validation design

The method was south validated for 55 pesticides, isomers or degradation products (appendix 1) in wheat, rye and oat. The validation was performed on 5-6 replicates at each of the three concentration levels. The concentration levels were 0.01, 0.02 and 0.1 mg/kg. A blank sample was included. The tests were done on different days by two different technicians.

4. Chromatograms and calibration curves

The calibration curve is determined by the analysis of each of the analysts at least 4 calibration levels, i.e. 0.003, 0.01, 0.033, 0.1 and 0.333 $\mu\text{g}/\text{ml}$. The calibration curves were best fitted to a linear curve. The quantification was performed from the mean of two calibration curves surrounding the samples. The majority of the correlation coefficients (R) were higher or equal to 0.98. Examples of chromatograms obtained when analysing the extracts by GC-MS/MS are presented in figure 1. Examples of calibration curves are presented in Figure 2.





Figure 1: Examples of chromatograms for cadusafos and ethoprophos obtained when analysing extract of wheat spiked with 0.1 mg/kg.

Compound name: Cadusafos
Correlation coefficient: $r = 0.999307$, $r^2 = 0.998614$
Calibration curve: $119913 * x + -69.5052$
Response type: External Std, Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



Compound name: Ethoprophos
Correlation coefficient: $r = 0.999023$, $r^2 = 0.998047$
Calibration curve: $78907 * x + -74.5637$
Response type: External Std, Area
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




Figure 2. Examples of calibration curves for cadusafos and ethoprophos (concentrations from 0.003-0.333 $\mu\text{g/ml}$)

5. Validation parameters

Precision – Repeatability

Repeatability was calculated for all pesticides and degradation products on all three spiking levels.

Repeatability is given as the relative standard deviation on the result from two or more analysis at the same sample, done by the same technician, on the same instrument and within a short period of time. Repeatability in this validation was calculated from the 5-6 replicate determinations.

Repeatability were calculated as given in ISO 5725-2².

Appendix 2 and 3 shows the relative repeatability for the validated pesticides, isomers and degradation products.

Accuracy – Recovery

The accuracy was determined by recovery, samples were spiked at three concentration levels. In appendix 2 and 3 recovery, repeatability and limit of quantification (LOQ) are given for the validated pesticides, isomers and degradation products for all three concentration levels (0.01 mg/kg, 0.02 mg/kg and 0.1 mg/kg). Recoveries may be seen in appendix 2 and 3.

Robustness

The QuEChERS method has earlier by Anastassiades et al. 2003¹ in connection with the development of the method been shown to be robust.

Limit of quantification, LOQ

Quantification limits (LOQ) are calculated from the results at the lowest accepted spike level, as 6 times the standard deviation (absolute recovery). The quantification limits are given in appendix 2 and 3.

6. Criteria for the acceptance of validation results

For the pesticides to be accepted as validated the following criteria for precision and trueness must to be fulfilled:

1. The relative standard deviation of the repeatability must be less than or equal to the standard deviation proposed by Horwitz³.
2. The average relative recovery must be between 70 and 120%⁴.

If the above mentioned criteria have been meet, the detection limits have been calculated.

7. Results and discussion

LC-MS-MS compounds

The multi-residue method has been tested for 50 pesticides, isomers and degradation products in wheat, rye, and oat using LC-MS/MS. The relative repeatability (RSD_r) varied between 2-22 %, however most of the values were around 6-14%. For the majority of the pesticides the recovery was in the range of 75-100% at all three concentration levels.

For wheat the criteria for acceptance were met for 32 out of 50 pesticides, isomers and degradation products (listed in Table 1). For rye and oat the criteria for acceptance were met for 34 out of 50 pesticides, isomers and degradation products (listed in Table 2-3). The lowest calibration level (LCL) was 0.0033 µg/ml and most of the LOQs are above 0.010 mg/kg.

The pesticides, isomers and degradation product which has been validated presented in table 1-3 are divided in to two groups, one for the compounds for which the acceptance criteria could be meet (Accepted) and those which could not meet the acceptance criteria (Not accepted).

The results for the different pesticides which were accepted for LC-MS/MS are listed in Appendix 2.

Table 1: Compounds validated and not accepted as validated for wheat.

Wheat		
Accepted (32 compounds)		
Avermectin B	Fenthionsulfon	Pencycuron
Boscalid	Fenthionsulfoxid	Phosmet
Bromuconazole sum	Fipronil	Phoxim
Cadusafos	Fipronil sulfone	Propagite
Clofentezin	Fipronil-sulfide	Spinosad A
Clothianidin	Formetanate	Tebufenozide
Ethoprophos	Fosthiazate	Teflubenzuron
Etofenprox	Lufenuron	Thiacloprid
Fenamiphos	Mepanipyrim	Thiamethoxam
Fenamiphos sulphone	Methoxyfenozide	Triazophos
Fenoxycarb	Paclobutrazole	
Not accepted (18 compounds)		
Amitraz	Fenthion	N-2,4-dimethylphenylformamide
Diafenthiurone	Flufenoxuron	2.4 dimethylphenyl-N-methylformamidine
Dinocap	Imidacloprid	Pyrethrins
Diphenylamine	Indoxacarb	Tolcloflos-methyl
Fenamidon	Prothioconazole	Zoxamide
Fenamiphos sulfoxide	Prothioconazole desthio	Spinosad D

Table 2: Compounds validated and not accepted as validated for oat.

Oat		
Accepted (34 compounds)		
Avermectin B	Fipronil	Propagite
Boscalid	Fipronil sulfone	Spinosad A
Bromuconazole sum	Fipronil-sulfide	Spinosad D
Cadusafos	Formetanate	Tebufenozide
Clofentezin	Fosthiazate	Teflubenzuron
Clothianidin	Lufenuron	Thiacloprid
Diphenylamine	Mepanipyrim	Thiamethoxam
Ethoprophos	Methoxyfenozide	Tolcloflos-methyl
Etofenprox	Paclobutrazole	Triazophos
Fenamiphos	Pencycuron	Zoxamide
Fenamiphos sulphone	Phosmet	
Fenoxycarb	Phoxim	
Not accepted (16 compounds)		
Amitraz	Fenthionsulfon	N-2,4-dimethylphenylformamide
Diafenthiurone	Fenthionsulfoxid	Prothioconazole desthio
Dinocap	Flufenoxuron	Pyrethrins
Fenamidon	Imidacloprid	2.4 dimethylphenyl-N-methylformamidine
Fenamiphos sulfoxide	Indoxacarb	
Fenthion	Prothioconazole	

Table 3: Compounds validated and not accepted as validated for rye.

Rye		
Accepted (34 compounds)		
Boscalid	Fipronil	Propagite
Bromuconazole sum	Fipronil sulfone	Spinosad A
Cadusafos	Fipronil-sulfide	Spinosad D
Clofentezin	Formetanate	Tebufenozide
Clothianidin	Fosthiazate	Teflubenzuron
Diphenylamine	Lufenuron	Thiacloprid
Ethoprophos	Mepanipyrim	Thiamethoxam
Etofenprox	Methoxyfenozide	Tolcloflos-methyl
Fenamiphos	Paclobutrazole	Triazophos
Fenamiphos sulphone	Pencycuron	Zoxamide
Fenoxycarb	Phosmet	
Fenthionsulfoxid	Phoxim	
Not accepted (16 compounds)		
Amitraz	Fenthion	N-2,4-dimethylphenylformamide
Avermectin B	Fenthionsulfon	2.4 dimethylphenyl-N-methylformamidine
Diafenthiurone	Flufenoxuron	Prothioconazole desthio
Dinocap	Imidacloprid	Pyrethrins
Fenamidon	Indoxacarb	
Fenamiphos sulfoxide	Prothioconazole	

GC-MS-MS compounds

The multi-residue method has been tested for 32 pesticides, isomers and degradation products in wheat, rye, and oat using GC-MS/MS. The relative repeatability (RSD_r) varied between 2-20 %, however most of the values were below 15 %. For the majority of the pesticides the recovery was in the range of 75-100% at all three concentration levels.

For wheat and oat the criteria for acceptance were met for 27 out of 32 pesticides, isomers and degradation products (listed in Table 4-5). For rye the criteria for acceptance were met for 28 out of 32 pesticides, isomers and degradation products (listed in Table 6).

The pesticides, isomers and degradation product which has been validated presented in table 4-6. The results for the different pesticides which were accepted are listed in Appendix 3.

Table 4: Compounds validated and not accepted as validated for wheat.

Wheat		
Accepted (27 compounds)		
Boscalid	Fenamiphos-sulfon	Indoxacarb
Bromuconazole sum	Fenthion	Mepanipyrim
Cadusafos	Fenthion-sulfon	Paclobutrazole
Chlorfenapyr	Fenthion-sulfoxid	Pencycuron
Clofentezin	Fipronil	Phosmet
Diphenylamine	Fipronil-sulfide	Tefluthrin
Ethoprofos	Fipronil-sulfone	Thiamethoxam
Etofenprox	Flufenoxuron	Tolclofos-methyl
Fenamiphos	Fosthiazate sum	Triazophos
Not accepted (5 compounds)		
Lufenuron	Propargite sum	Amitraz
Formetanate	Fenamiphos-sulfoxide	

Table 5: Compounds validated and not accepted as validated for oat.

Oat		
Accepted (27 compounds)		
Boscalid	Fenamiphos-sulfoxide	Lufenuron
Bromuconazole sum	Fenthion	Mepanipyrim
Cadusafos	Fenthion-sulfon	Paclobutrazole
Chlorfenapyr	Fenthion-sulfoxid	Pencycuron
Clofentezin	Fipronil	Phosmet
Ethoprofos	Fipronil-sulfide	Tefluthrin
Etofenprox	Fipronil-sulfone	Thiamethoxam
Fenamiphos	Flufenoxuron	Tolclofos-methyl
Fenamiphos-sulfon	Indoxacarb	Triazophos
Not accepted (5 compounds)		
Amitraz	Formetanate	Propargite sum
Diphenylamine	Fosthiazate sum	

Table 6: Compounds validated and not accepted as validated for rye.

Rye		
Accepted (28 compounds)		
Boscalid	Fenamiphos-sulfoxide	Mepanipyrim
Bromuconazole sum	Fenthion	Paclobutrazole
Cadusafos	Fenthion-sulfon	Pencycuron
Chlorfenapyr	Fenthion-sulfoxid	Phosmet
Clofentezin	Fipronil	Tefluthrin
Diphenylamine	Fipronil-sulfide	Thiamethoxam
Ethoprophos	Fipronil-sulfone	Tolclofos-methyl
Etofenprox	Flufenoxuron	Triazophos
Fenamiphos	Fosthiazate sum	
Fenamiphos-sulfon	Indoxacarb	
Not accepted (4 compounds)		
Amitraz	Lufenuron	Propargite sum
Formetanate		

8. Conclusions

In conclusion 44 of 55 pesticides, isomers and degradation products were validated for the QuEChERS method using GC-MS/MS and LC-MS/MS for the analysis. 25 compounds were validated on both GC-MS/MS and LC-MS/MS. 13 compounds were only validated on LC-MS/MS and 6 compounds were only validated on GC-MS/MS.

9. References

- 1 <http://www.quechers.com/> or Anastassiades et al., J. AOAC Int., vol. 86, no. 2, p. 412, 2003
- 2 ISO 5725-2:1994. Accuracy (trueness and precision) of measurement methods and results – Part 2. Basic method for the determination of repeatability and reproducibility of standard measurement method. First edition. December 1994.
- 3 W. Horwitz, Anal. Chem., 1982; 54, 67A.
- 4 Method Validation and Quality Control Procedures for Pesticide Residue Analysis in Food and Feed, Document No SANCO/10684/2010, 01/01/2010, European Commission, Brussels, 2010.

Appendix 1. MRM transitions for the pesticides sought validated.

GC-MS/MS		Parent 1	Fragment 1	Col. energi 1	Intensity	Parent 2	Fragment 2	Col. energi 2	Intensity
1	Clofentezin	137	102	10	53600	102	75	10	58800
2	Lufenuron	352.9	203	15	677000	173.9	110	10	329000
3	Diphenylamine	169	168	10	7750000	168	167	10	6450000
4	Ethoprophos	158	97	15	12300000	200	158	5	13400000
5	Pencycuron	180	125	10	7880000	125	89	10	1370000
6	Cadusafos	159	131	5	18800000	158	114	5	6680000
7	Flufenoxuron	331	268	15	2440000	268	241	15	923000
8	Tefluthrin	177.1	127	15	18800000	197.1	141	10	6350000
9	Tolcloflos-methyl	265	250	15	14300000	267	252	10	5550000
10	Fenthion	278	109	15	10100000	278	124	20	3870000
11	Thiamethoxam	212	139	10	3170000	246.9	182	10	1680000
12	Fosthiazate 1	195	103	5	1460000	283	103	15	303000
13	Fosthiazate 2	195	103	5	1460000	283	103	15	303000
14	Fipronil-sulfide	419.9	351	10	8720000	350.9	255	15	12100000
15	Fipronil	367	213	20	2950000	367	255	15	1110000
16	Formatanate	163	122	10	415000	122	77	10	124000
17	Paclobutrazole	236	125	15	1430000	238	127	15	564000
18	Mepanipirim	222.1	207	15	3590000	223.1	207	20	1880000
19	Fenamiphos	154	139	10	1310000	303	195	5	2000000
20	Fipronil sulfone	382.9	255	20	2270000	255	228	10	660000
21	Chlorfenapyr	327.9	247	15	536000	363.9	248	20	142000
22	Fenthionsulfoxid	277.9	109	15	1480000	246.9	125	10	856000
23	Fenthionsulfon	309.9	105	10	4010000	136.1	91	20	481000
24	Triazophos	257	162	5	11200000	285	162	10	10600000
25	Fenamiphos sulfoxide	304	122	15	3380000	288.0	243.0	5.0	1880000
26	Fenamiphos sulphone	320	292	10	11700000	291.9	134	15	3350000
27	Phosmet	160	77	20	1820000	160	133	10	1300000
28	Amitraz	162.2	132	10	3090000	132.1	117	10	3130000

GC-MS/MS		Parent 1	Fragment 1	Col. energi 1	Intensity	Parent 2	Fragment 2	Col. energi 2	Intensity
29	Boscalid	341.9	140	15	6840000	167	139	20	295000
30	Etofenprox	163	107	15	14700000	107	77	15	1340000
31	Indoxacarb	264	176	10	922000	527	264	5	427000
32	Propagite 1	135	107	5	1380000	350.2	81	10	300000
33	Propagite 2	135	107	5	1380000	350.2	81	10	300000
34	Bromuconazole 1	295	173	10	1260000	293	173	15	1010000
35	Bromuconazole 2	295	173	10	1260000	293	173	15	1010000

LC-MS/MS ESI-		Moder ion	CV	Daughter-1	CE	Daughter-2	CE
1	Clothianidin	247.923	50	58.134	13	150	13
2	Methoxyfenozide	367.43	42	149.001	23	105.114	31
3	Tebufenozide	351.43	42	149.035	23	105.044	31
4	Fipronil	435.22	42	330.152	13	250.069	23
5	Fipronil-sulfide	419.215	42	261.999	13	383.161	23
6	Fipronil sulfone	451.18	42	282.066	31	244.064	40
7	Prothioconazole	342.176	20	99.932	23	124.996	23
8	Lufenuron	509.216	13	339.186	13	175.009	31
9	Diafenthurone	383.916	50	191.432	30	268.677	25
10	Teflubenzuron	379.176	42	339.107	23	195.953	13
11	Dinocap	295.1	20	209.1	20	192.1	20

LCMSMS ESI+		Moder ion	CV	Daughter-1	CE	Daughter-2	CE
1	Formetanate	222.226	30	165.04	20	92.92	20
2	Thiamethoxam	292.256	50	211.21	13	181.169	23
3	2,4 dimethylphenyl-N-methylformamidine	163.715	10	122.955	20	132.936	20
4	Imidacloprid	256.41	30	175.27	13	209.39	13
5	Thiacloprid	253.204	10	126.115	23	90.287	40
6	N-2,4-dimethylphenylformamide	150.25	34	106.827	35	132.213	20
7	Fenamiphos sulfoxide	320.365	50	108.133	13	292.323	40
8	Fenamiphos sulphone	353.441	50	188.224	31	266.266	23
9	Fenthionsulfoxid	295.1	25	108.9	20	279.7	15
10	Fenthionsulfon	311.1	25	125	20	108.9	20
11	Fosthiazate	284.355	30	104.093	25	228.145	13

LCMSMS ESI+		Moder ion	CV	Daughter-1	CE	Daughter-2	CE
12	Phosmet	335.314	20	161.202	31	160	39
13	Fenamidon	312.484	30	236.329	40	92.117	13
14	Boscalid	343.35	57	307.22	20	140.25	20
15	Paclobutrazole	294.331	20	70.23	31	125.153	31
16	Triazophos	314.312	42	162.236	13	119.193	40
17	Ethoprophos	243.241	50	97.049	23	131.056	31
18	Diphenylamine	170	25	93.1	20	151.9	20
19	Mepanipyrim	224.204	42	77.196	23	106.17	39
20	Fenamiphos	304.41	50	217.079	23	202.069	31
21	Prothioconazole desthio	312.484	50	70.34	35	125.424	30
22	Fenoxycarb	302.331	42	88.215	31	116.228	13
23	Fenthion	279.1	25	169.1	20	247.1	20
24	Bromuconazole 1	378.232	42	159.072	23	70.182	23
25	Bromuconazole 2	378.232	42	159.072	23	70.182	23
26	Zoxamide	336.39	36	159.125	20	187.19	40
27	Phoxim	299.298	10	115.029	31	77.1	23
28	Indoxacarb	528.458	42	203.494	31	150.34	31
29	Tolcloflos-methyl	301	25	175	20	268.9	20
30	Clofentezin	303.1	20	138.1	20	102.1	20
31	Pencycuron	329.458	50	125.123	23	218	13
32	Spinosad A	733.01	30	142.182	20	98	20
33	Cadusafos	271.331	50	159.025	13	97.018	40
34	Spinosad D	746.5	30	142.311	40	98.168	45
35	Flufenoxuron	489.614	22	158.147	13	140.949	40
36	Propagite	368.478	30	175.231	13	231.406	13
37	Pyrethrins	329.503	34	128.211	31	162.138	13
38	Amitraz	294.486	30	163.181	13	122.104	31
39	Avermectin B	890.5	41	145.2	43	305.1	35
40	Etofenprox	394.2	25	177.3	20	107.1	20

Appendix 2. Repeatability, recovery and limit of quantification for compounds validated on LC-MS-MS.

In the tables are presented repeatability and LOQ for compounds validated on LC-MS-MS. Values outside the acceptance criteria are marked by grey. The matrixes are wheat, oat and rye.

WHEAT	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	
Clothianidin	70	26	103	3	95	3	0.012
Methoxyfenozide	59	30	87	11	91	11	0.011
Tebufenozide	64	8	91	8	101	7	0.009
Fipronil	69	11	87	11	91	8	0.009
Fipronil-sulfide	80	28	96	12	98	9	0.015
Fipronil sulfone	65	36	96	6	102	8	0.007
Lufenuron	57	13	92	9	86	6	0.009
Teflubenzuron	53	30	79	11	90	10	0.010
Formetanate	88	21	84	19	87	3	0.012
Thiamethoxam	127	57	109	36	86	8	0.041
Thiacloprid	93	13	92	9	86	6	0.008
Avermectin B	99	62	113	63	96	22	0.133
Etofenprox	70	13	69	20	78	12	0.006
Fenamiphos sulphone	87	18	88	26	15	15	0.011
Fenthionsulfoxid	98	44	105	3	100	7	0.029
Fenthionsulfon	104	27	124	40	97	21	0.126
Fosthiazate	84	21	85	9	73	10	0.012
Phosmet	88	12	94	14	96	10	0.007
Boscalid	102	35	75	26	94	15	0.023
Paclobutrazole	83	31	98	5	106	8	0.017
Triazophos	80	20	80	12	77	12	0.010
Ethoprophos	79	19	102	9	99	9	0.010
Mepanipyrim	114	64	77	108	102	19	0.118
Fenamiphos	85	17	82	17	74	8	0.010
Fenoxycarb	106	4	98	21	79	7	0.003
Bromuconazole sum	112	22	90	10	76	15	0.016
Phoxim	95	18	82	23	95	18	0.011
Propagite	89	14	98	19	90	11	0.008
Clofentezin	52	99	122	79	76	10	0.047
Pencycuron	87	27	81	17	80	10	0.015
Spinosad A	81	23	100	16	81	10	0.013
Cadusafos	91	15	97	6	85	10	0.009

OAT	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	
Clothianidin	89	15	109	17	95	17	0.009
Methoxyfenozide	88	16	98	18	103	18	0.009
Tebufenozide	69	20	97	20	97	6	0.023
Fipronil	81	18	98	18	103	5	0.015
Fipronil-sulfide	85	29	91	16	91	9	0.016
Fipronil sulfone	87	21	94	14	91	6	0.012
Lufenuron	68	14	120	13	102	7	0.018
Teflubenzuron	71	18	77	16	78	3	0.008
Formetanate	84	20	92	10	83	9	0.011
Thiamethoxam	96	30	111	17	82	19	0.019
Thiacloprid	95	14	120	13	102	7	0.009
Avermectin B	60	29	83	18	94	22	0.018
Etofenprox	71	16	86	29	93	13	0.007
Fenamiphos sulphone	93	14	111	5	111	6	0.008
Fosthiazate	92	11	118	5	108	8	0.007
Phosmet	93	16	107	4	105	6	0.010
Boscalid	85	20	108	9	108	6	0.011
Paclobutrazole	93	20	117	8	108	7	0.012
Triazophos	95	12	99	19	104	6	0.007
Ethoprophos	88	16	99	11	101	8	0.009
Diphenylamine	53	51	59	51	104	6	0.040
Mepanipyrim	100	14	99	17	101	8	0.009
Fenamiphos	87	11	103	7	100	8	0.006
Fenoxycarb	79	16	104	6	106	6	0.008
Bromuconazole sum	82	10	107	10	107	5	0.005
Zoxamide	947	204	97	22	97	11	0.025
Phoxim	2092	179	89	22	96	10	0.023
Propagite	86	15	99	15	106	7	0.009
Tolcloflos-methyl	972	206	87	34	97	18	0.111
Clofentezin	273	146	107	14	100	20	0.019
Pencycuron	97	9	108	4	105	6	0.005
Spinosad A	75	16	87	19	102	9	0.008
Cadusafos	84	10	103	7	100	10	0.006
Spinosad D	83	32	105	15	105	8	0.018

RYE	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	
Clothianidin	89	18	109	7	97	11	0.011
Methoxyfenozide	77	16	97	12	90	13	0.008
Tebufenozide	73	34	92	10	93	9	0.011
Fipronil	87	21	97	12	90	13	0.012
Fipronil-sulfide	77	19	107	8	102	9	0.010
Fipronil sulfone	84	12	103	4	96	12	0.007
Lufenuron	81	17	80	5	81	11	0.009
Teflubenzuron	112	17	76	14	76	13	0.012
Formetanate	89	15	110	8	103	9	0.009
Thiamethoxam	55	36	82	18	106	9	0.018
Thiacloprid	81	18	101	9	104	9	0.010
Etofenprox	70	34	79	16	80	9	0.015
Fenamiphos sulphone	82	25	105	10	110	6	0.014
Fenthionsulfoxid	63	37	88	19	106	10	0.019
Fosthiazate	77	11	106	12	111	7	0.005
Phosmet	80	20	102	5	104	6	0.010
Boscalid	91	38	102	10	110	8	0.012
Paclobutrazole	69	18	110	6	110	9	0.008
Triazophos	62	23	100	10	90	13	0.012
Ethoprophos	80	18	112	14	111	8	0.009
Diphenylamine	57	116	48	19	90	13	0.07
Mepanipyrim	73	26	80	16	82	19	0.013
Fenamiphos	76	22	103	6	108	5	0.011
Fenoxycarb	91	20	100	11	110	9	0.012
Bromuconazole sum	91	15	110	14	108	7	0.009
Zoxamide	60	84	121	60	117	21	0.152
Phoxim	76	53	100	23	112	5	0.027
Propagite	71	12	91	20	85	8	0.005
Tolcloflos-methyl	64	170	97	86	93	4	0.023
Clofentezin	109	50	106	14	108	7	0.018
Pencycuron	78	8	93	12	92	8	0.004
Spinosad A	70	22	79	13	89	5	0.010
Cadusafos	70	19	98	7	106	4	0.009
Spinosad D	52	56	67	10	84	10	0.054

Appendix 3. Repeatability, recovery and limit of quantification for compounds validated on GC-MS-MS.

WHEAT	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	LOQ
Clofentezin	73	10	71	7	75	7	0.005
Diphenylamine	96	11	98	11	70	12	0.004
Ethoprophos	102	5	103	5	76	10	0.012
Pencycuron	92	9	88	7	76	19	0.005
Cadusafos	93	10	105	3	86	9	0.006
Flufenoxuron	73	16	80	13	67	11	0.008
Tefluthrin	93	16	103	8	109	6	0.011
Tolclofos-methyl	98	11	102	6	108	9	0.007
Fenthion	96	14	99	4	98	8	0.009
Thiamethoxam	97	12	103	8	101	9	0.008
Fosthiazate sum	101	3	111	5	111	2	0.002
Fipronil-sulfide	93	18	103	5	103	8	0.011
Fipronil	102	16	103	8	109	6	0.011
Paclobutrazole	100	16	110	5	102	9	0.010
Mepanipyrim	95	7	100	7	105	12	0.004
Fenamiphos	89	10	94	8	82	19	0.006
Fipronil-sulfone	92	11	101	5	110	4	0.007
Chlorfenapyr	97	26	91	13	107	5	0.017
Fenthion-sulfoxid	85	18	92	6	107	10	0.010
Fenthion-sulfon	84	18	96	4	108	4	0.010
Triazophos	94	11	100	7	113	1	0.007
Fenamiphos-sulfon	97	15	96	8	100	5	0.010
Bromuconazole sum	85	17	100	9	105	12	0.009
Phosmet	102	9	90	9	74	14	0.006
Boscalid	87	14	91	6	101	7	0.008
Etofenprox	85	0	81	6	73	19	0.005
Indoxacarb	85	8	101	10	106	6	0.005

OAT	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	LOQ
Clofentezin	105	15	89	9	102	9	0.011
Lufenuron	101	5	71	4	70	4	0.003
Ethoprophos	115	3	120	3	81	7	0.010
Pencycuron	78	10	77	6	82	6	0.005
Cadusafos	119	12	119	2	82	6	0.010
Flufenoxuron	86	15	82	8	70	6	0.008
Tefluthrin	103	15	110	4	81	2	0.012
Tolclofos-methyl	111	14	119	5	76	5	0.010
Fenthion	85	16	88	5	70	7	0.009
Thiamethoxam	112	8	116	1	76	3	0.006

Fipronil-sulfide	118	14	119	8	88	3	0.011
Fipronil	126	15	131	4	81	2	0.011
Paclobutrazole	124	20	147	2	89	4	0.016
Mepanipyrim	111	31	116	7	70	5	0.023
Fenamiphos	99	12	110	3	74	3	0.008
Fipronil-sulfone	114	13	119	1	91	3	0.010
Chlorfenapyr	102	13	102	18	80	3	0.009
Fenthion-sulfoxid	132	18	144	4	84	13	0.068
Fenthion-sulfon	100	13	110	2	76	7	0.009
Triazophos	101	11	110	4	75	8	0.007
Fenamiphos-sulfoxide	105	9	106	3	72	8	0.006
Fenamiphos-sulfon	97	21	95	3	74	15	0.013
Bromuconazole sum	99	6	103	3	81	2	0.004
Phosmet	92	11	80	9	70	4	0.007
Boscalid	99	10	110	4	80	2	0.007
Etofenprox	95	0	102	5	73	4	0.010
Indoxacarb	91	11	92	2	77	17	0.006

RYE	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	Spike level mg/kg	Horwitz, %	LOQ
	0.011	32	0.02	29	0.104	22	
	Recovery, %	RSDr, %	Recovery, %	RSDr, %	Recovery, %	RSDr, %	
Clofentezin	92	18	103	9	99	9	0.011
Diphenylamine	107	7	110	7	81	3	0.009
Ethoprophos	98	6	118	6	89	4	0.008
Pencycuron	98	10	102	7	82	2	0.007
Cadusafos	96	12	117	2	84	2	0.007
Flufenoxuron	76	18	85	7	109	6	0.009
Tefluthrin	91	8	110	4	83	4	0.005
Tolclofos-methyl	95	12	114	3	84	3	0.007
Fenthion	96	11	119	3	79	3	0.007
Thiamethoxam	85	15	100	8	79	7	0.008
Fosthiazate sum	99	12	108	7	97	31	0.008
Fipronil-sulfide	98	13	115	2	83	3	0.008
Fipronil	101	8	110	4	83	4	0.005
Paclobutrazole	100	11	119	5	80	6	0.007
Mepanipyrim	99	11	110	5	75	9	0.007
Fenamiphos	91	11	110	8	76	4	0.007
Fipronil-sulfone	113	22	119	2	99	5	0.016
Chlorfenapyr	95	14	116	3	81	2	0.009
Fenthion-sulfoxid	94	13	88	3	77	12	0.008
Fenthion-sulfon	97	13	102	4	83	4	0.008
Triazophos	87	16	109	2	78	9	0.009
Fenamiphos-sulfoxide	100	7	89	10	87	16	0.005
Fenamiphos-sulfon	95	11	94	5	87	22	0.007
Bromuconazole sum	98	12	110	3	78	2	0.007
Phosmet	96	20	90	8	83	5	0.013
Boscalid	93	13	109	6	80	4	0.008
Etofenprox	90	0	107	8	74	2	0.008
Indoxacarb	87	11	95	8	275	63	0.006