

Matrix effects impact using different chromatographic column sizes

Table of contents

| | |
|--|----|
| 1. Aim and scope | 1 |
| 2. Short description | 1 |
| 3. Apparatus and consumables | 1 |
| 4. Chemicals | 1 |
| 5. Procedure | 2 |
| 6. Results | 3 |
| 7. References | 9 |
| APPENDIX I: Mass Transitions and retention times | 10 |

1. Aim and scope

This document reports the matrix effects evaluation data in two chromatographic columns for 209 pesticides included in the European Union Multi Annual Control Program (EU-MACP) and the Working Document [1,2] using a multiresidue method by LC-MS/MS. These columns, which differ in their particle size and length, were two of the most popular column sizes used in routine laboratories.

2. Short description

The retention times for 209 pesticides were experimentally determined for both chromatographic columns tested. Then, homogenous tomato, orange, and onion samples were extracted using QuEChERS without a clean-up step. The obtained extracts were used to prepare matrix-matched calibration curves with concentrations ranging between 5 µg/L and 500 µg/L. The four calibration curves (including solvent) were injected in the LC-MSMS instrument using each chromatographic column and matrix effects were evaluated.

3. Apparatus and consumables

- Automatic pipettes, suitable for handling volumes from 10 µL to 5000 µL and from 1 mL to 5 mL.
- Graduated 10 mL pipette.
- 50 mL and 15 mL PTFE centrifuge tubes.
- Vortex Shaker IKATM 4 Basic.
- Axial shaker Agytax SR1 CP57.
- Centrifuge Orto Alresa Consul 21, suitable for the centrifuge tubes employed in the procedure and capable of achieving at least 4000 rpm.
- Concentration workstation.
- Injection vials, 2 mL, suitable for LC and GC auto-sampler.

4. Chemicals

- Acetonitrile ultra-gradient grade
- Trisodium citrate dihydrate
- Disodium hydrogenocitrate sesquihydrate
- Sodium chloride
- Anhydrous magnesium sulphate
- Anhydrous calcium chloride
- Ammonium formate

- Ultra-pure water
- Methanol HPLC grade
- Formic acid
- Pesticide standards

5. Procedure

5.1. Sample preparation

Tomato, orange and onion samples were extracted following the QuEChERS extraction procedure without a clean-up step. Removing the clean-up step the amount of matrix components in the extract increase.

5.2. Pesticide stock solutions and working mix solutions

Individual pesticide stock solutions (1000–2000 mg/L) were prepared in acetonitrile or ethyl acetate and were stored in screw-capped glass vials in the dark at -20 °C. Working mixes were prepared in 10 mL volumetric flasks by pipetting the appropriate volume of each stock solution.

5.3. Instrumentation and analytical conditions for the LC- MS/MS system

5.3.1. Nexera UC (Shimadzu)

- Mobile phase A: Water (0.1 % formic acid, 5 mM ammonium formate, 2 % MeOH)
- Mobile phase B: Methanol (0.1 % formic acid, 5 mM ammonium formate, 2 % water)
- Column temperature: 40 °C
- Flow rate: Column A: 0.3 mL/min
Column B: 0.5 mL/min
- Injection volume: 2 µL

These parameters were used in combination with the two different chromatographic columns evaluated in this technical report:

- Column A: fully porous C8 2.1x100 mm and 1.8 µm particle size
- Column B: fully porous C8 2.1x150 mm and 3.0 µm particle size

Same run time was applied in both columns evaluation (15min)
 Mobile phase gradient:

| min | Mobile phase A | Mobile phase B |
|-----|----------------|----------------|
| 0 | 100 % | 0 % |
| 1.5 | 75 % | 25 % |
| 10 | 0 % | 100 % |
| 12 | 0 % | 100 % |

Re-equilibration time with initial mobile phase set for 3 minutes.

5.3.2. 8060 triple quadrupole system (Shimadzu)

- Ionisation mode ESI: Positive mode and negative mode
- Capillary (positive and negative): 4 kV
- Nebulizing gas flow: 3 L/min
- Heating gas flow: 10 L/min
- Drying gas flow: 10 L/min
- Interface temperature: 350°C
- DL temperature: 250°C
- Heat block temperature: 400°C

6. Results

6.1. Matrix effects evaluation

Matrix effects were calculated comparing the slopes of the matrix-matched calibration curves (5-500 µg/Kg range) with the slope of the calibration curve prepared in solvent. The equation used to calculate this was as follows:

$$ME (\%) = \left(\left(\frac{\text{matrix - matched slope}}{\text{solvent slope}} \right) - 1 \right) \times 100$$

Coeluting compounds from the matrix can produce enhancement or suppression of the analyte signal. When the analytes enter the ion source at the same time as the matrix compounds, usually signal suppression is the predominant phenomenon in liquid chromatography. Matrix effects modify the reproducibility and accuracy of the method leading to an erroneous quantification. Matrix effects can differ

depending on the column size used. The elution of matrix compounds could not be the same if different particle size/length columns are used. Between 0 and 20% is considered a low or non-existent matrix effect; however, modifications of the signal between 20 and 50% and >50% are considered medium and strong matrix effects, respectively.

Table 1. Matrix effects in three different matrices using both columns

| Matrix Effects | | Tomato | Orange | Onion |
|----------------|-------|--------|--------|-------|
| Column A | 0-20 | 90% | 70% | 52% |
| | 20-50 | 9% | 20% | 43% |
| | >50 | 1% | 10% | 5% |
| Column B | 0-20 | 88% | 61% | 55% |
| | 20-50 | 5% | 24% | 34% |
| | >50 | 7% | 15% | 11% |




Figure 1. Matrix effects in tomato matrix.




Figure 2. Matrix effects in orange matrix.




Figure 3. Matrix effects in onion matrix.

6.2. Conclusions

The objective of this study was to compare two of the typical columns employed in routine laboratories and evaluate the matrix effects by analysing three different matrices (tomato, orange, and onion). These matrix effects can be affected by the distribution of the coextracted compounds along with the chromatogram and the number of plates provided by the column. The column dimensions were 2.1x100mm, 1.8 μm particle size (Column A) and 2.1x150mm, 3.5 μm particle size (Column B). Both columns used fully porous particles and employed C8 as stationary phase. Matrix effects results were similar in both columns. However, results were slightly better using column A, specially considering the percentage of compounds with strong matrix effects (above 50%). This fact makes sense taking into account chromatography resolution. Narrower peaks might be able to avoid the coelution with matrix components reducing ion suppression in the source (Figure 4). Other concerns should take into account like solvent waste which increase by 60% using column B as higher mobile phase flow was necessary. Nonetheless, despite in general terms matrix effects were similar, sometimes one compound can produce a different response in a specific matrix depending on the column, as can be observed in figures 5-8 with chlorantraniliprole and bromacil in orange matrix.




Figure 4. Chromatograms of coumaphos at the concentration level of 10 $\mu\text{g}/\text{Kg}$ in tomato matrix.




Figure 5. Solvent and matrix-matched slopes for chlorantraniliprole using column A.




Figure 6. Solvent and matrix-matched slopes for chlorantraniliprole using column B.




Figure 7. Solvent and matrix-matched slopes for Bromacil using column A.




Figure 8. Solvent and matrix-matched slopes for Bromacil using column B.

7. References

[1] COMMISSION IMPLEMENTING REGULATION (EU) 2019/533 of 28 March 2019 concerning a coordinated multiannual control programme of the Union for 2020, 2021 and 2022 to ensure compliance with maximum residue levels of pesticides and to assess the consumer exposure to pesticide residues in and on food of plant and animal origin.

[2] Working document on pesticides to be considered for inclusion in the national control programmes to ensure compliance with maximum residue levels of pesticides residues in and on food of plant and animal origin (SANCO/12745/2013).

APPENDIX I: MASS TRANSITIONS
Table 1. Detection parameters for the selected compounds analysed by LC-MS/MS.

| Compound | Precursor Ion (1) | Product Ion (1) | Precursor Ion (2) | Product Ion (2) | Ret. Time Column A (min) | Ret. Time Column B (min) |
|---------------------|-------------------|-----------------|-------------------|-----------------|--------------------------|--------------------------|
| 2,4-D | 220.9 | 163.05 | 218.9 | 161.15 | 7.762 | 8.815 |
| Acephate | 184 | 143 | 184 | 125 | 3.341 | 4.265 |
| Acetamiprid | 223.1 | 126.1 | 223.1 | 99.05 | 5.524 | 6.749 |
| Alachlor | 270.1 | 162.15 | 270.1 | 147.1 | 9.359 | 10.423 |
| Ametoctradin | 276.1 | 176 | 276.1 | 149.15 | 10.028 | 10.89 |
| Anilofos | 368 | 125 | 368 | 199.05 | 9.71 | 10.645 |
| Atrazine | 216.1 | 174.1 | 216.1 | 104.05 | 7.968 | 9.388 |
| Azinphos-ethyl | 346 | 97 | 346 | 132.1 | 9.101 | 10.429 |
| Azinphos-methyl | 318 | 76.95 | 318 | 132.05 | 8.266 | 9.744 |
| Azoxystrobin | 404 | 371.95 | 404 | 328.95 | 8.491 | 9.795 |
| Bendiocarb | 224.1 | 109 | 224.1 | 81.05 | 7.248 | 8.559 |
| Bifenazate | 301.1 | 198.1 | 301.1 | 170 | 9.063 | 10.175 |
| Bifenazate_diazene | 299 | 213.05 | 299 | 184 | 9.916 | 11.23 |
| Bitertanol | 338 | 99.15 | 338 | 269.05 | 9.807 | 10.743 |
| Boscalid | 343 | 306.95 | 343 | 271.95 | 8.718 | 9.791 |
| Bromacil | 261 | 204.9 | 261 | 188 | 7.212 | 8.443 |
| Bromuconazole | 375.9 | 159 | 379.9 | 161.05 | 9.004 | 9.995 |
| Bupirimate | 317.2 | 166 | 317.2 | 108 | 9.283 | 10.443 |
| Buprofezin | 306.2 | 201.05 | 306.2 | 116.15 | 10.539 | 11.405 |
| Carbaryl | 202.1 | 127 | 202.1 | 145.1 | 7.436 | 8.821 |
| Carbendazim | 192.1 | 160.15 | 192.1 | 132.15 | 4.359 | 5.858 |
| Carbofuran | 222.1 | 123.15 | 222.1 | 165 | 7.274 | 8.51 |
| Chlorantraniliprole | 483.9 | 452.9 | 483.9 | 285.9 | 8.292 | 9.411 |
| Chlorbromuron | 295 | 205.9 | 295 | 182.1 | 8.703 | 10 |
| Chlorfenvinphos | 359 | 169.95 | 359 | 155.15 | 9.812 | 10.692 |
| Chlorfluazuron | 539.9 | 382.85 | 539.9 | 158 | 10.91 | 11.721 |
| Chloridazon | 222.1 | 104.1 | 222.1 | 92.15 | 5.511 | 7.009 |
| Chlorotoluron | 213.1 | 72.15 | 213.1 | 139.95 | 7.758 | 9.076 |
| Chloroxuron | 291.1 | 72.15 | 291.1 | 218.05 | 8.979 | 10.108 |
| Chlorpyrifos | 350 | 197.95 | 352 | 199.85 | 10.578 | 11.651 |
| Chromafenozide | 394.9 | 175.15 | 394.9 | 91.1 | 9.268 | 10.26 |
| Clofentezine | 303 | 138.15 | 303 | 102.1 | 9.659 | 11.135 |
| Clomazone | 239.9 | 125 | 241.9 | 127 | 8.448 | 9.723 |
| Coumaphos | 363 | 227 | 363 | 306.8 | 9.576 | 10.805 |
| Cyazofamid | 325 | 108.1 | 325 | 261 | 9.282 | 10.377 |
| Cyflufenamid | 413.1 | 295.05 | 413.1 | 241 | 9.974 | 10.934 |
| Cymoxanil | 199.1 | 111.15 | 199.1 | 128.15 | 5.867 | 7.129 |
| Cyproconazole | 292.1 | 125.05 | 292.1 | 89 | 9.163 | 10.061 |

| | | | | | | |
|--------------------------|-------|--------|-------|--------|--------|--------|
| Cyprodinil | 226.1 | 93 | 226.1 | 108 | 9.265 | 10.73 |
| Cyromazine | 167.2 | 85.1 | 167.2 | 125.15 | 2.101 | 3.789 |
| Deet | 192.1 | 91 | 192.1 | 119 | 8.15 | 9.307 |
| Demeton-S-methyl | 231.1 | 89 | 231.1 | 61.1 | 7.303 | 8.533 |
| Demeton-S-methyl-sulfone | 263 | 169.05 | 263 | 109.05 | 4.462 | 5.532 |
| Diazinon | 305.1 | 169.1 | 305.1 | 96.95 | 9.792 | 10.9 |
| Dichlorvos | 238 | 220.9 | 238 | 109.1 | 7.184 | 8.473 |
| Dicrotophos | 237.9 | 127 | 237.9 | 112.15 | 4.951 | 6.06 |
| Diethofencarb | 268.2 | 226.05 | 268.2 | 124.15 | 8.517 | 9.69 |
| Difenoconazole | 406.1 | 250.9 | 406.1 | 188 | 9.948 | 10.87 |
| Difenoxyuron | 287.1 | 72.05 | 287.1 | 123.15 | 7.956 | 9.36 |
| Diflubenzuron | 311 | 158.1 | 311 | 141.05 | 9.282 | 10.392 |
| Dimethoate | 230 | 125 | 230 | 198.9 | 5.595 | 6.909 |
| Dimethomorph | 388.1 | 301 | 388.1 | 165.1 | 8.815 | 9.84 |
| Diniconazole | 326.1 | 70.05 | 328 | 70 | 10.07 | 10.902 |
| Diuron | 233 | 72.1 | 233 | 46.15 | 7.762 | 9.043 |
| Desethyl_terbuthylazine | 202 | 146.05 | 202 | 104 | 7.482 | 8.73 |
| Edifenphos | 311 | 109 | 311 | 111.05 | 9.631 | 10.765 |
| Emamectin B1a | 886.4 | 158.2 | 886.4 | 82.05 | 10.231 | 10.789 |
| EPN | 324.1 | 295.85 | 324.1 | 156.95 | 9.996 | 11.174 |
| Epoconazole | 330 | 121.1 | 330 | 101.1 | 9.209 | 10.257 |
| Ethiofencarb | 226.1 | 107 | 226.1 | 77 | 7.688 | 8.934 |
| Ethion | 385 | 143 | 385 | 198.9 | 10.579 | 11.562 |
| Ethiprole | 397 | 350.9 | 397 | 254.85 | 8.719 | 9.587 |
| Ethirimol | 210.2 | 140.2 | 210.2 | 98.15 | 6.564 | 7.692 |
| Ethoprophos | 243.1 | 172.9 | 243.1 | 131 | 9.353 | 10.365 |
| Etofenprox | 394.2 | 177.05 | 394.2 | 107.05 | 11.443 | 12.331 |
| Etoxazole | 360.1 | 141.1 | 360.1 | 113.05 | 10.809 | 11.725 |
| Fenamidone | 312.1 | 92.1 | 312.1 | 236 | 8.756 | 9.919 |
| Fenamiphos | 304.1 | 216.95 | 304.1 | 201.95 | 9.445 | 10.375 |
| Fenamiphos-sulfone | 336.1 | 266 | 336.1 | 188 | 7.385 | 8.35 |
| Fenamiphos-sulfoxide | 320.1 | 233 | 320.1 | 108 | 7.257 | 8.185 |
| Fenarimol | 331 | 268 | 331 | 189 | 9.223 | 10.286 |
| Fenazaquin | 307.2 | 161.1 | 307.2 | 131.2 | 11.031 | 11.964 |
| Fenbuconazole | 337.1 | 70.1 | 337.1 | 125.05 | 9.299 | 10.218 |
| Fenhexamid | 302.1 | 97.1 | 302.1 | 143 | 9.162 | 10.087 |
| Fenobucarb | 208.1 | 95 | 208.1 | 152.1 | 8.66 | 9.744 |
| Fenoxycarb | 302.1 | 116.15 | 302.1 | 256.05 | 9.393 | 10.518 |
| Fenpropimorph | 304.2 | 117 | 304.2 | 147.1 | 8.496 | 9.35 |
| Fenpyroximate | 422.1 | 366 | 422.1 | 138.1 | 10.876 | 11.798 |
| Fenthion-sulfone | 311 | 125.1 | 311 | 109 | 7.544 | 8.898 |

| | | | | | | |
|----------------------|-------|--------|-------|--------|--------|--------|
| Fenthion-sulfoxide | 295 | 279.9 | 295 | 109 | 7.373 | 8.692 |
| Fenuron | 165 | 72.15 | 165 | 77 | 5.378 | 6.855 |
| Fipronil | 435 | 330 | 435 | 250.05 | 9.522 | 10.336 |
| Flazasulfuron | 408 | 182.05 | 408 | 139 | 8.26 | 9.441 |
| Fluacrypyrim | 427.2 | 145.15 | 427.2 | 205.05 | 10.121 | 11.105 |
| Fluazifop | 327.7 | 282 | 327.7 | 254 | 8.63 | 9.593 |
| Fludioxonil | 247 | 180.15 | 247 | 169.15 | 8.747 | 10.063 |
| Flufenacet | 364.1 | 152.05 | 364.1 | 194 | 9.325 | 10.262 |
| Flufenoxuron | 489 | 140.9 | 489 | 158.1 | 10.729 | 11.52 |
| Fluopicolide | 383 | 172.85 | 383 | 144.95 | 8.939 | 10.021 |
| Fluopyram | 397 | 145 | 397 | 208 | 9.217 | 10.161 |
| Fluquinconazole | 376 | 307 | 376 | 349 | 9.041 | 10.132 |
| Flusilazole | 316.1 | 247 | 316.1 | 165.1 | 9.451 | 10.356 |
| Flutriafol | 302.1 | 70.05 | 302.1 | 123 | 7.97 | 9.06 |
| Fluxapyroxad | 382 | 342 | 382 | 313.9 | 8.885 | 9.917 |
| Fosthiazate | 284.1 | 104.1 | 284.1 | 227.85 | 7.772 | 8.896 |
| Haloxfop | 362.1 | 315.95 | 364.1 | 318.15 | 9.468 | 10.321 |
| Hexaconazole | 314.1 | 70 | 314.1 | 159 | 9.896 | 10.7 |
| Hexythiazox | 353.1 | 228 | 353.1 | 168 | 10.718 | 11.554 |
| Imazalil | 297.1 | 159.05 | 297.1 | 200.85 | 7.698 | 8.618 |
| Imidacloprid | 256.1 | 174.95 | 256.1 | 209 | 5.02 | 6.261 |
| Indoxacarb | 527.7 | 203 | 527.7 | 218 | 10.098 | 10.966 |
| Ioxynil | 369.8 | 127.05 | 369.8 | 215 | 8.08 | 9.31 |
| Iprodione | 330.1 | 245 | 332.1 | 247 | 9.351 | 10.262 |
| Iprovalicarb | 321.2 | 119.15 | 321.2 | 203 | 9.314 | 10.242 |
| Isoprocarb | 194.1 | 95 | 194.1 | 137 | 8.041 | 9.225 |
| Isoprothiolane | 290.8 | 188.9 | 290.8 | 231 | 8.875 | 10.168 |
| Isoproturon | 207.2 | 72.15 | 207.2 | 165 | 8.069 | 9.257 |
| Isoxaflutole | 360.1 | 251 | 360.1 | 144 | 8.158 | 9.337 |
| Kresoxim-methyl | 314.1 | 235 | 314.1 | 267 | 9.55 | 10.758 |
| Lenacil | 234.9 | 153.15 | 234.9 | 136 | 8.042 | 9.234 |
| Linuron | 248.8 | 182.05 | 248.8 | 160 | 8.57 | 9.881 |
| Malathion | 348.1 | 127.05 | 348.1 | 330.9 | 8.939 | 10.126 |
| Malathion D10 | 340.9 | 132.1 | 340.9 | 290.05 | 8.902 | 10.111 |
| Mandipropamid | 412.1 | 327.9 | 412.1 | 125.1 | 8.763 | 9.841 |
| Mepanipyrim | 224.1 | 106.05 | 224.1 | 104.1 | 8.947 | 10.507 |
| Metaflumizone | 506.8 | 178.05 | 506.8 | 287 | 10.444 | 11 |
| Metalaxyl | 280.1 | 220 | 280.1 | 192.05 | 8.147 | 9.232 |
| Metconazole | 320.1 | 70.15 | 320.1 | 125.05 | 9.837 | 10.648 |
| Methamidophos | 142.2 | 93.95 | 142.2 | 124.9 | 2.662 | 3.92 |
| Methidathion | 320 | 145 | 320 | 302.8 | 8.239 | 9.683 |
| Methiocarb | 226.1 | 121.1 | 226.1 | 169.05 | 8.681 | 9.885 |
| Methiocarb-sulfoxide | 242.1 | 185.05 | 242.1 | 122.1 | 5.309 | 6.623 |
| Methomyl | 163 | 87.9 | 163 | 106.15 | 4.306 | 5.753 |

| | | | | | | |
|----------------------|-------|--------|-------|--------|--------|--------|
| Methoxyfenozide | 369.2 | 149.15 | 369.2 | 91.15 | 9.08 | 10.127 |
| Metobromuron | 259 | 170 | 259 | 148.1 | 7.802 | 9.197 |
| Metolachlor | 284.1 | 252.05 | 284.1 | 176.1 | 9.422 | 10.472 |
| Metolcarb | 166.1 | 109.1 | 166.1 | 94.05 | 6.807 | 8.203 |
| Metrafenone | 409 | 209.1 | 409 | 227 | 9.917 | 11.033 |
| Monocrotophos | 240.9 | 193 | 240.9 | 127.1 | 4.658 | 5.791 |
| Monolinuron | 215.1 | 99.1 | 215.1 | 148 | 7.558 | 8.989 |
| Monuron | 199.1 | 72.15 | 199.1 | 126 | 6.933 | 8.328 |
| Myclobutanil | 289.1 | 70.05 | 289.1 | 125 | 9.04 | 9.923 |
| Neburon | 274.8 | 87.95 | 274.8 | 114.15 | 9.548 | 10.576 |
| Nitenpyram | 271.1 | 225 | 271.1 | 237 | 4.112 | 5.009 |
| Novaluron | 493 | 141.05 | 493 | 158 | 10.18 | 11.074 |
| Omethoate | 214.1 | 125 | 214.1 | 183 | 3.641 | 4.574 |
| Oxadialgyl | 358.1 | 340.9 | 358.1 | 223 | 9.864 | 10.858 |
| Oxadixyl | 296.2 | 279.05 | 296.2 | 219.05 | 6.696 | 7.85 |
| Oxamyl | 237.1 | 72.1 | 237.1 | 90 | 4.109 | 5.2 |
| Oxasulfuron | 407.1 | 150.15 | 407.1 | 107 | 6.801 | 8.156 |
| Paclobutrazol | 294.1 | 70.1 | 294.1 | 125.05 | 8.924 | 9.86 |
| Penconazole | 284.1 | 158.95 | 284.1 | 123 | 9.67 | 10.549 |
| Pencycuron | 329.1 | 125 | 329.1 | 89 | 10.016 | 10.933 |
| Pendimethalin | 282.2 | 212 | 282.2 | 194 | 10.648 | 11.713 |
| Phenthoate | 321 | 79.05 | 321 | 246.85 | 9.548 | 10.78 |
| Phosalone | 368 | 182 | 368 | 111 | 9.841 | 10.848 |
| Phosmet | 335 | 160.1 | 335 | 317.9 | 8.325 | 9.773 |
| Phoxim | 299 | 77.1 | 299 | 129.1 | 9.799 | 11.069 |
| Pirimicarb | 239.2 | 182.05 | 239.2 | 85.05 | 6.863 | 8.427 |
| Pirimicarb-desmethyl | 225.1 | 168.05 | 225.1 | 180 | 5.188 | 6.5 |
| Pirimiphos-methyl | 306.1 | 164.1 | 306.1 | 108.15 | 9.852 | 11.069 |
| Prochloraz | 376 | 307.95 | 376 | 70 | 9.72 | 10.66 |
| Profenofos | 375 | 304.7 | 375 | 346.75 | 10.343 | 11.255 |
| Promecarb | 208.1 | 109.1 | 208.1 | 151.15 | 8.896 | 10.005 |
| Prometryn | 242.1 | 158 | 242.1 | 200.15 | 8.885 | 10.205 |
| Propamocarb | 189.2 | 102.15 | 189.2 | 74.15 | 3.724 | 4.598 |
| Propaquizafop | 443.8 | 100.15 | 443.8 | 299 | 10.325 | 11.241 |
| Propargite | 368.2 | 231.1 | 368.2 | 175.1 | 10.811 | 11.613 |
| Propazine | 230.2 | 146.15 | 230.2 | 188.1 | 8.677 | 9.942 |
| Propiconazole | 342 | 158.9 | 342 | 69.1 | 9.776 | 10.628 |
| Propoxur | 209.9 | 111.15 | 209.9 | 93.1 | 7.21 | 8.448 |
| Propyzamide | 256 | 190 | 256 | 173 | 8.937 | 10.096 |
| Proquinazid | 373.1 | 288.8 | 373.1 | 330.9 | 10.831 | 11.889 |
| Prosulfocarb | 252.1 | 91 | 252.1 | 128.1 | 10.278 | 11.273 |
| Prothiophos | 345 | 240.9 | 347 | 242.75 | 11.157 | 12.087 |
| Pymetrozine | 218.1 | 105 | 218.1 | 78.1 | 3.418 | 4.291 |
| Pyraclostrobin | 388 | 194.1 | 388 | 133 | 9.693 | 10.884 |

| | | | | | | |
|----------------------------------|-------|--------|-------|--------|--------|--------|
| Pyridaben | 365.2 | 147.1 | 365.2 | 308.9 | 11.116 | 11.876 |
| Pyridaphenthion | 341.1 | 188.95 | 341.1 | 204.9 | 9.013 | 10.161 |
| Pyridate | 379.1 | 206.95 | 379.1 | 104.15 | 11.283 | 12.022 |
| Pyrimethanil | 200.1 | 107 | 200.1 | 183 | 8.251 | 9.832 |
| Pyriproxyfen | 322.1 | 96.05 | 322.1 | 184.95 | 10.452 | 11.613 |
| Quinalphos | 299.1 | 163 | 299.1 | 147.05 | 9.453 | 10.805 |
| Quinoclamine | 208 | 88.95 | 208 | 105.05 | 6.673 | 8.362 |
| Quinoxyfen | 308 | 197 | 308 | 271.95 | 10.52 | 11.536 |
| Quizalofop (free acid) | 345 | 299 | 345 | 243.95 | 9.184 | 10.11 |
| Quizalofop-ethyl | 373.1 | 298.9 | 373.1 | 270.9 | 10.196 | 11.176 |
| Rotenone | 395.1 | 213 | 395.1 | 192 | 9.27 | 10.519 |
| Simazine | 202 | 104.05 | 202 | 96.05 | 7.13 | 8.713 |
| SpinetoramA | 748.5 | 142.15 | 748.5 | 98.05 | 10.012 | 10.55 |
| SpinetoramB | 760.6 | 142.1 | 760.6 | 98.15 | 10.227 | 10.793 |
| Spinosad A | 732.6 | 142.2 | 732.6 | 98.1 | 9.674 | 10.283 |
| Spinosad D | 746.6 | 142.1 | 746.6 | 98.1 | 9.934 | 10.538 |
| Spirodiclofen | 411.1 | 313.05 | 411.1 | 294.9 | 11.051 | 11.775 |
| Spiromesifen | 371.2 | 273.1 | 371.2 | 255.2 | 10.905 | 11.509 |
| Spirotetramat | 374.1 | 216 | 374.1 | 302 | 9.179 | 10.2 |
| Tebuconazole | 308.2 | 70.05 | 308.2 | 125.05 | 9.653 | 10.479 |
| Tebufenozide | 353.2 | 133.1 | 353.2 | 297 | 9.606 | 10.48 |
| Tebufenpyrad | 334.2 | 117 | 334.2 | 145.15 | 10.436 | 11.237 |
| Teflubenzuron | 378.8 | 339 | 378.8 | 358.9 | 10.393 | 11.357 |
| Terbutylazine | 230.1 | 174.05 | 230.1 | 104 | 8.787 | 10.052 |
| Terbutylazine-desethyl | 202.1 | 146.05 | 202.1 | 104.05 | 7.482 | 8.728 |
| Terbutryn | 242.1 | 157.95 | 242.1 | 90.9 | 8.887 | 10.205 |
| Tetraconazole | 372 | 159 | 372 | 70.1 | 9.293 | 10.181 |
| Tetramethrin | 332.2 | 164.1 | 332.2 | 135.1 | 10.47 | 11.374 |
| Thiabendazole | 201.8 | 175 | 201.8 | 131.15 | 4.804 | 6.658 |
| Thiacloprid | 252.8 | 126.05 | 252.8 | 90.1 | 5.964 | 7.244 |
| Thiamethoxam | 292 | 211.1 | 292 | 132 | 4.449 | 5.562 |
| Thiobencarb | 257.8 | 125.1 | 257.8 | 89 | 9.935 | 10.988 |
| Triazophos | 314.1 | 162.05 | 314.1 | 119.15 | 9.025 | 10.293 |
| Trichlorfon | 257 | 109 | 257 | 220.8 | 5.598 | 6.742 |
| Tricyclazole | 190.1 | 136 | 190.1 | 109 | 6.17 | 7.699 |
| Trifloxystrobin | 408.8 | 186.1 | 408.8 | 145.1 | 10.169 | 11.195 |
| Triflumizole | 346 | 278 | 348 | 280 | 10.267 | 11.099 |
| Triflumuron | 359 | 156.05 | 359 | 139.05 | 9.79 | 10.769 |
| Triticonazole | 318.1 | 70.15 | 318.1 | 125 | 9.213 | 10.07 |
| Tritosulfuron | 446 | 194.9 | 446 | 145.05 | 8.341 | 9.39 |
| XMC (3,5-xyllyl methylcarbamate) | 180.1 | 123.1 | 180.1 | 108 | 7.671 | 8.795 |
| Zoxamide | 336 | 186.95 | 336 | 159 | 9.756 | 10.731 |