

Shimadzu solutions for helium shortage in GC/GC-MS

Nerea Lorenzo Parodi
Product Specialist GC-MS
Shimadzu Europa GmbH

Index

- Introduction
- Alternative carrier gases
- Safety features with Hydrogen
- Application examples with Hydrogen
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Helium crisis

- Global demand is increasing
- Maintenance work in existing plants
- Delayed opening of new plant



- Significant price increase
- Unavailability

Index


- Introduction
- **Alternative carrier gases**
- Safety features with Hydrogen
- Application examples with Hydrogen
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Gases in gas chromatography

Carrier gases **Pros (+)** / **Cons (-)**

	Helium	Nitrogen	Hydrogen
Diffusivity	+	-	++



Gases in gas chromatography

Carrier gases **Pros (+)** / **Cons (-)**

	Helium	Nitrogen	Hydrogen
Diffusivity	+	-	++
Viscosity	-	-	+
Cost	--	+	+
Availability	--	+	+
Inertness	+	+	-
Others	For some detection techniques (e.g. BID) it is the only solution	Peak resolution worsens with increasing velocities	Cannot be used with some detectors, explosion risk (!)
Conclusion	Ideal carrier gas	Alternative to helium for easy to separate components	For many applications, good results and faster

Index


- Introduction
- Alternative carrier gases
- **Safety features with Hydrogen**
- Application examples with Hydrogen
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Safe GC operation with Hydrogen

Safety features of Shimadzu GCs

1. Dangerous leaks are automatically detected by **fast responding AFC**
2. Smaller leaks affecting GC result are detected by **automated leak check**
3. Additional safety by **AFC/APC hydrogen option** – reduces maximum hydrogen flow/pressure – and **optional hydrogen sensor** – continuous monitoring of hydrogen content in oven air – .




GC-2030 built-in hydrogen sensor:

Hydrogen concentration in GC oven:

1 % \Rightarrow Error Message, AFC flow stop

2 % \Rightarrow System shuts down

Explosion level of hydrogen is 4% in air



Index

- Introduction
- Alternative carrier gases
- Safety features with Hydrogen
- **Application examples with Hydrogen**
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Application examples with Hydrogen

- Pesticide Residue Analysis with GC-MS/MS
 - 120 pesticides, 5 - 500 ppb
 - Apple extract, QuEChERs for calibration
- 40 % time saved with good resolution




Figure 2: Fluvalinate-1 and -2 signals for helium (left) and hydrogen (right)

Application examples with Hydrogen




- Enough sensitivity at 5 ppb, Intensity drop compound dependent

Methacrifos




Sulfotep

alpha-HCH

He



H₂



Index


- Introduction
- Alternative carrier gases
- Safety features with Hydrogen
- Application examples with Hydrogen
- **Reducing He consumption**
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Cost reduction with Gas Saver and Ecology Mode


Gas Saver:

- In the method file, after a set time, the split ratio can be reduced to save carrier gas



Ecology Mode:

- A shutdown method file can be set with, e.g. lower temperatures and flows
- Several temperatures can be turned off
- Ecology mode can be activated:
 - Manually at any time
 - Automatically at the end of a batch



Cost reduction with Gas Saver and Ecology Mode

Savings example:

Normal Method

Carrier Gas: He	
Flow Control Mode:	Linear Velocity
Pressure:	69.4 kPa
Total Flow:	31.5 mL/min
Column Flow:	1.22 mL/min
Linear Velocity:	40.0 cm/s
Purge Flow:	6.0 mL/min
Split Ratio:	20.0

- Total Flow = 31.5 mL / min
- 50 L He Bottle (200 bar)
- Total “run” time = 220 days

Ecology Method

Carrier Gas: He	
Flow Control Mode:	Linear Velocity
Pressure:	32.7 kPa
Total Flow:	6.5 mL/min
Column Flow:	0.75 mL/min
Linear Velocity:	31.4 cm/s
Purge Flow:	2.0 mL/min
Split Ratio:	5.0



- Total Flow = 6.5 mL / min
- 50 L He Bottle (200 bar)
- Total “run” time = 1068 days

If used 10 h a day, 7 days a week the “lifetime” of the He bottle **doubles!**

Cost reduction with Gas Selector

Gas Selector:

- The gas can be switched with a **few clicks** in the software
- Either gas can be used for the analysis
- After the analysis, the GC can **automatically switch** to the alternative carrier gas
- The system can **automatically shutdown/re-start**
- Due to efficient purge function, **10 minutes** are enough for GC to restart with full helium separation performance



Automated shutdown and restart function of GC-2030

Index

- Introduction
- Alternative carrier gases
- Safety features with Hydrogen
- Application examples with Hydrogen
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- **Conclusion**
- Resources



Conclusion


- Shimadzu offers several measures for the Helium Crisis:

- Use of alternative carrier gases
- Be safe while using Hydrogen thanks to:

- Fast responding AFC
- Automatic leak check

- Reduce Helium consumption thanks to:

- Gas saver
- Ecology Mode
- Gas selector



- Hydrogen Sensor
- AFC/APC Hydrogen option



Index

- Introduction
- Alternative carrier gases
- Safety features with Hydrogen
- Application examples with Hydrogen
- Reducing He consumption
 - Gas Saver and Ecology Mode
 - Gas selector
- Conclusion
- Resources



Resources

[Carrier Gas Theory](#)

[Short on helium? Consider the carrier-gas alternatives](#)

[Countermeasures and Solutions for Helium Gas Supply Shortages](#)

[Alternative GCMS Carrier Gas](#)

[Hydrogen as carrier gas](#)

[Gas Chromatograph Hydrogen Gas Safety](#)

[Nexis GC-2030](#)



Thank you for your attention!!

Questions?

info@shimadzu.eu



SHIMADZU

Excellence in Science