

## Impact of the different calibration approaches in the quality of the results

Table of contents

<b>1. Aim and scope</b>	<b>2</b>
<b>2. Short description</b>	<b>2</b>
<b>3. Apparatus and consumables</b>	<b>2</b>
<b>4. Chemicals</b>	<b>2</b>
<b>5. Procedure</b>	<b>3</b>
<b>6. Results</b>	<b>6</b>
<b>7. Conclusions</b>	<b>28</b>
<b>8. References</b>	<b>29</b>
<b>APPENDIX I: MASS TRANSITIONS</b>	<b>30</b>

## 1. Aim and scope

This document reports the impact that the different ways of **calibration** for quantification, described in document SANTE/11312/2021 [1], have in the quality of the results for **171** pesticides included in the European Union Multi Annual Control Program (EU-MACP) and the Working Document [2,3] using a multiresidue method by LC-MS/MS and GC-MS/MS.

## 2. Short description

**Different ways** of calibration have been tested to assess the influence they have in the **quality** of the results: multiple- and single-level matrix matched calibration using the same matrix, multiple-level calibration of a different matrix of the same group, procedural calibration, sample standard addition and extract standard addition. The different calibration approaches were tested in one representative commodity of the three main FV commodity groups: **high water** content, **high acid** and water content, and **high oil** and very low water content. The extracts were prepared using **QuEChERS** as the extraction method.

## 3. Apparatus and consumables

- Automatic pipettes, suitable for handling volumes from 1 µL to 5000 µL and from 1 mL to 5 mL.
- Graduated 10 mL pipette.
- 50 mL and 15 mL PTFE centrifuge tubes.
- Vortex Shaker IKATM 4 Basic.
- Axial shaker Agytax SR1 CP57.
- Centrifuge Orto Alresa Consul 21, suitable for the centrifuge tubes employed in the procedure and capable of achieving at least 4000 rpm.
- Concentration workstation.
- Injection vials, 2 mL, suitable for LC and GC auto-sampler.

## 4. Chemicals

- Acetonitrile ultra-gradient grade
- Trisodium citrate dihydrate
- Disodium hydrogenocitrate sesquihydrate
- Sodium chloride
- Anhydrous magnesium sulphate
- Anhydrous calcium chloride
- Primary secondary amine (PSA)
- Supel QuE Z-Sep

- Ammonium formate
- Ultra-pure water
- Methanol HPLC grade
- Formic acid
- Pesticide standards

## 5. Procedure

### 5.1. Sample preparation

The main matrices used in the study were **tomato**, **orange** and **avocado**. Additionally, multiple-level calibration in **potato**, **strawberry** and **olive extract** were prepared.

### 5.2. Pesticide stock solutions and working mix solutions

Individual pesticide stock solutions (1000–2000 mg/L) were prepared in acetonitrile or ethyl acetate and were stored in screw-capped glass vials in the dark at -20 °C. Working mixes were prepared in 10 mL volumetric flasks by pipetting the appropriate volume of each stock solution.

### 5.3. Extraction procedure

#### QuEChERS

1. Weigh 10 g ± 0.1 g of sample in 50 mL PTFE centrifuge tube.
2. Add 10 mL of acetonitrile and 10 µL of 10 mg/L carbendazim-d3, malathion-d10 and TPP (procedure internal standards).
3. Shake the sample using an automatic axial shaker for 4 min.
4. Add 4 g of magnesium sulphate, 1 g of sodium chloride, 1 g of trisodium citrate dihydrate and 0.5 g of disodium hydrogenocitrate sesquihydrate.
5. Shake the samples again in the automatic shaker for 4 min.
6. Centrifuge the tubes at 3700 rpm for 5 min.
7. Transfer 5 mL of the supernatant to a 15 mL PTFE tube containing:
  - a. 750 mg magnesium sulphate and 125 mg PSA for matrices with high water content and high acid content and high water content.
  - b. 750 mg magnesium sulphate and 175 mg Z-Sep for matrices with high oil content and very low water content.
8. Vortex the tube for 30 sec.
9. Centrifuge the tubes at 3700 rpm for 5 min.
10. Add 40 µL of formic acid 5% in acetonitrile to option a in step 7.
11. Analysis:

- a. LC-MS/MS: dilute 100 mL extract with 400 mL of water containing dimethoate-d6 at 0.050 mg/L (Injection Internal Standard).
- b. GC-MS/MS: recovery samples were evaporated under a gentle N<sub>2</sub> current and reconstituted with ethyl acetate prior to injection.  
Standard calibration levels were prepared using the blank extract and reconstitute with the pesticide mix in ethyl acetate.  
In both cases, add 0.002 mL of Lindane-d6 (Injection Internal Standard) and transfer the contents to a 0.20 mL spring-loaded insert.

#### 5.4. Measurement

LC system was operated in multiple reaction monitoring mode (MRM). Selected reaction monitoring (SRM) experiments were carried out to obtain the maximum sensitivity for the detection of the target molecules. For confirmation of the studied compounds, two SRM transitions and a correct ratio between the abundances of the two optimized SRM transitions (SRM2/SRM1) were used, along with retention time matching. The mass transitions used are presented in Appendix I.

#### 5.5. Instrumentation and analytical conditions for the LC- MS/MS system

##### 5.5.1. 1290 UHPLC (Agilent)

- Column: Zorbax Eclipse Plus C8 2.1x100 mm and 1.8 µm particle size (Agilent)
- Mobile phase A: Water (0.1% formic acid, 5mM ammonium formate, 2% MeOH)
- Mobile phase B: Methanol (0.1% formic acid, 5mM ammonium formate, 2% H<sub>2</sub>O)
- Column temperature: 35°C
- Flow rate: 0.3 mL/min
- Injection volume: 5 µL.

Mobile phase gradient for pesticides analysed

min	Mobile phase A	Mobile phase B
0	100 %	0 %
1.5	75 %	25 %
10	0 %	100 %
12	0 %	100 %

Re-equilibration with initial phase: 2.5 minutes

#### 5.5.2. 6490 triple quadrupole system (Agilent)

- Ionisation mode: Positive mode and negative mode
- Capillary (positive and negative): 3000 V
- Nebulizer: 45 psi
- Nozzle: 400 V
- Drying gas flow: 13 L/min
- Drying gas temperature: 120°C
- Sheath gas flow: 10 L/min
- Sheath gas temperature: 375°C
- High Pressure RF (positive): 150 V
- High Pressure RF (negative): 110 V
- Low Pressure RF (positive): 60 V
- Low Pressure RF (negative): 60 V

### 5.6. Instrumentation and analytical conditions for the GC- MS/MS system

#### 5.6.1. Intuvo 9000 GC system (Agilent)

- Columns: 2 planar columns HP-5MS UI (15 m long × 0.25 mm i.d. × 0.25 µm film thickness)
- Injection mode: Splitless
- Sample injection volume: 1 µL
- Inlet temperature: 80 °C hold for 0.1 min, then up to 300 °C at 600 °C/min, hold for 5 min and then to 250 °C at 100 °C/min
- Carrier gas: Helium at constant flow = 1.28 mL/min column 1, 1.48 mL/min column 2
- Carrier gas purity: 99.999 %
- Oven temperature: 60 °C for 0.5 min, up to 170 °C at 80 °C/min, and up to 310 °C at 20 °C/min (hold for 3.5 min)
- Post Run: 2.1 min, 310 °C

#### 5.6.2. 7410 triple quadrupole system (Agilent)

- Ionisation mode: electron impact ionisation
- Temperature of the transfer line: 280 °C
- Temperature of ion source: 280 °C
- Collision gas: nitrogen
- Collision gas purity: 99.999 %
- Solvent delay: 2.6 minutes

## 6. Results

In order to evaluate the impact that different ways of quantification had in the quality of the results, 10 µg/kg and 100 µg/kg analytical standard mixes were prepared in blank extracts of tomato, orange and avocado, as representative commodities from the three studied groups. Those were considered as the "fake samples" that should be quantified using the different quantification strategies. This way, each one of those fake samples was quantified as follows:

- Sample of **tomato** (10 and 100 µg/kg) quantified with a **calibration curve in tomato extract**.
- Sample of **tomato** (10 and 100 µg/kg) quantified with a **single-point calibration** using each one of the levels of the calibration curve (5 µg/kg, 10 µg/kg, 50 µg/kg, 100 µg/kg and 200 µg/kg) in **tomato extract**.
- Sample of **tomato** (10 and 100 µg/kg) quantified with a **calibration curve in potato extract**.
- Recovery test (10 and 100 µg/kg) in **tomato** quantified by **procedural calibration in tomato**.
- Sample of **tomato** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the extract.
- Recovery test of **tomato** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the sample.

- Sample of **orange** (10 and 100 µg/kg) quantified with a **calibration curve in orange extract**.
- Sample of **orange** (10 and 100 µg/kg) quantified with a **single-point calibration** using each one of the levels of the calibration curve (5 µg/kg, 10 µg/kg, 50 µg/kg, 100 µg/kg and 200 µg/kg) in **orange extract**.
- Sample of **orange** (10 and 100 µg/kg) quantified with a **calibration curve in strawberry extract**.
- Recovery test (10 and 100 µg/kg) in **orange** quantified by **procedural calibration in orange**.
- Sample of **orange** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the extract.
- Recovery test of **orange** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the sample.

- Sample of **avocado** (10 and 100 µg/kg) quantified with a **calibration curve in avocado extract**.
- Sample of **avocado** (10 and 100 µg/kg) quantified with a **single-point calibration** using each one of the levels of the calibration curve (5 µg/kg, 10 µg/kg, 50 µg/kg, 100 µg/kg and 200 µg/kg) in **avocado extract**.

- Sample of **avocado** (10 and 100 µg/kg) quantified with a **calibration curve** in **olive** extract.
- Recovery test (10 and 100 µg/kg) in **avocado** quantified by **procedural calibration in avocado**.
- Sample of **avocado** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the extract.
- Recovery test of **avocado** (10 µg/kg) quantified by **standard addition** (10 µg/kg, 20 µg/kg, 30 µg/kg and 40 µg/kg) to the sample.

Percentage variation of concentrations from the theoretical concentrations at the 10 µg/kg level

After quantification of all the samples by the methods described above, the **percentage variation** of the concentrations with respect to the theoretical ones (10 and 100 µg/kg) was calculated using the following equation:

$$\text{Variation (\%)} = \left| \left( \frac{\text{Concentration obtained } (\mu\text{g/kg})}{\text{Theoretical concentration } (\mu\text{g/kg})} \right) \times 100 - 100 \right|$$

The **results** obtained for the quantification of the samples at the 10 µg/kg level are shown in **Table 1**, whereas the results for the quantification of the 100 µg/kg level are presented in **Table 2**.

Pesticides with a variation of | 0-21% | are shown in **green**, with variations of | 21-25% | in **orange** and for cases with variations higher than | 25% | in **red**.

**Table 3** shows the percentage of analysed compounds that fall into each of the classifications at 10 µg/kg (green, orange or red). Those tested at 100 µg/kg are displayed in **Table 4**.

**Table 1.** Percentage variation of concentrations from the theoretical concentrations at the 10 µg/kg level. Pesticides with a variation of | 0-21% | are shown in green, | 21-25% | in orange and > | 25% | in red.

**A** = multiple-level calibration of the same matrix (tomato, orange and avocado); single-level calibration: **B** = 5 µg/kg, **C** = 10 µg/kg, **D** = 50 µg/kg, **E** = 100 µg/kg, **F** = 200 µg/kg; **G** = multiple-level calibration of another matrix from the same group of matrices (potato, strawberry and olive); **H** = procedural calibration; **I** = sample standard addition; **J** = extract standard addition

Pesticide	10 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J
2,4'-DDE										Red																			Red	
2-Phenylphenol										Yellow									Red	Red									Green	
Acephate																				Red									Red	
Acetamiprid																													Red	
Acrinathrin										Yellow																				
Aldicarb											Red																			
Aldicarb-sulfone																	Red				Red	Red								
Ametoctradin										Yellow																				Yellow
Azinphos-methyl																	Red				Yellow									Red
Azoxystrobin																			Yellow	Red										
Benalaxyd										Red										Red										
Benzovindiflupyr																		Red	Red											
Bifenthrin																	Red	Red												
Bitertanol										Yellow	Yellow	Red							Red											
Boscalid																			Yellow	Red										
Bromopropylate										Red																		Yellow	Red	

Pesticide	10 µg/kg																														
	TOMATO										ORANGE										AVOCADO										
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	
Bupirimate																															
Buprofezin																															
Carbaryl																															
Carbendazim																															
Carbofuran																															
Chlorantraniliprole																															
Chlorfenapyr																															
Chlorfluazuron																															
Chloridazon																															
Chlorpropham																															
Chlorpyrifos																															
Chlorpyrifos-methyl																															
Clomazone																															
Clopyralid																															
Cyflufenamid																															
Cyfluthrin																															
Cyhalofop-butyl																															
Cymoxanil																															
Cypermethrin																															
Cyproconazole																															
Cyprodinil																															
Cyromazine																															

Pesticide	10 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J
Deltamethrin										Y																				
Diazinon																					Y									
Dichloran																														
Dichlorvos																														
Dicofol, o, p'-																														
Dicofol, p, p'-										Y																				
Dieldrin																														
Diethofencarb																														
Difenoconazole																														
Diflubenzuron																														
Dimethoate																														
Dimethomorph																														
Diniconazole																														
Dinotefuran																														
Diuron																														
Emamectin B1a benzoate										Y																				
Emamectin B1b benzoate																														
Endosulfan sulfate																														
Endosulfan-alpha																														
Endosulfan-beta																														
Epoxiconazole																														

Pesticide	10 µg/kg																														
	TOMATO										ORANGE										AVOCADO										
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	
Ethion																															
Ethirimol																															
Ethoprophos																															
Etofenprox																															
Etoxazole																															
Famoxadone																															
Fenamidone																															
Fenamiphos																															
Fenamiphos-sulfone																															
Fenamiphos-sulfoxide																															
Fenarimol																															
Fenazaquin																															
Fenbuconazole																															
Fenhexamid																															
Fenitrothion																															
Fenobucarb																															
Fenoxy carb																															
Fenpicoxamid																															
Fenpropathrin																															
Fenpropidin																															
Fenpropimorph																															

Pesticide	10 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J
Fenpyrazamine																														
Fenpyroximate																														
Fenthion																														
Fenthion-sulfone																														
Fenthion-sulfoxide																														
Fenvalerate																														
Fipronil																														
Flonicamid																														
Fluazifop-p-butyl																														
Flubendiamide																														
Fludioxonil																														
Fluensulfone																														
Flufenacet																														
Flufenoxuron																														
Fluopicolide																														
Fluopyram																														
Fluquinconazole																														
Flusilazole																														
Flutriafol																														
Fluxapyroxad																														
Formetanate Hydrochloride																														
Fosthiazate																														

Pesticide	10 µg/kg																																						
	TOMATO										ORANGE										AVOCADO																		
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J									
Haloxyfop	Y	R				Y	R	R	R		Y				R	R	R																						
Haloxyfop-methyl															R		Y				Y	R	R	R	R	R													
Heptachlor				Y	R																																		
Hexaconazole							R												R																				
Hexythiazox															R	R	R																						
Imazalil																																							
Indoxacarb												R																											
loxynil																	R																						
Iprovalicarb																			Y																				
Isocarbophos						Y						R																											
Isoprothiolane																		R	Y																				
Isopyrazam												R																											
Isoxaflutole																		R																					
Kresoxim-methyl																				Y																			
Lambda-Cyhalothrin																			R																				
Linuron																																							
Lufenuron																		R	Y																				
Malathion																				Y																			
Mandipropamid																																							
Mepanipyrim																																							
Metaflumizone	R	R	Y	Y																																			
Metalaxyl																																							

Pesticide	10 µg/kg																														
	TOMATO										ORANGE										AVOCADO										
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	
Metconazole																															
Methamidophos																															
Methidathion																															
Methiocarb																															
Methomyl																															
Methoxyfenozide																															
Metrafenone																															
Molinate																															
Monocrotophos																															
Myclobutanil																															
Novaluron																															
Omethoate																															
Oxadiargyl																															
Oxadixyl																															
Oxamyl																															
Oxsulfuron																															
Oxathiapipronil																															
Paclobutrazol																															
Parathion-methyl																															
Penconazole																															
Pencycuron																															
Pendimethalin																															
Penflufen																															

Pesticide	10 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J
Penthiopyrad																														
Permethrin																														
Phenthroate																														
Phosmet																														
Picolinafen																														
Pirimicarb																														
Pirimicarb-desmethyl																														
Pirimiphos-methyl																														
Propaquizafop																														
Prothifos																														
Pyriofenone																														
Quinalphos																														
Quinoclamine																														
Quintozene																														
Quizalofop																														
Quizalofop-ethyl																														
Rotenone																														
Tetramethrin																														
Tolfenpyrad																														
Triallate																														
Triflumizole																														
Triticonazole																														

Pesticide	10 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J	A	B	C	D	E	F	G	H	I	J
Tritosulfuron	green	green	green	green	green	green	green	green	yellow		green	green	green	green	red	yellow	green	green	red	red	green	red	red	red	red	green	green	red	red	red
Zoxamide	red	green	green	green	red	yellow	green	green	green		green	green	green	green	green	red	green	green	yellow	green	green	green	green	green	green	green	green	yellow		

**A** = multiple-level calibration of the same matrix (tomato, orange and avocado); single-level calibration: **B** = 5 µg/kg, **C** = 10 µg/kg, **D** = 50 µg/kg, **E** = 100 µg/kg, **F** = 200 µg/kg; **G** = multiple-level calibration of another matrix from the same group of matrices (potato, strawberry and olive); **H** = procedural calibration; **I** = sample standard addition; **J** = extract standard addition

**Table 2.** Percentage variation of concentrations from the theoretical concentrations at the 100 µg/kg level Pesticides with a variation of | 0-21% | are shown in green, | 21-25% | in orange and > | 25% | in red.

**A** = multiple-level calibration of the same matrix (tomato, orange and avocado); single-level calibration: **B** = 5 µg/kg, **C** = 10 µg/kg, **D** = 50 µg/kg, **E** = 100 µg/kg, **F** = 200 µg/kg; **G** = multiple-level calibration of another matrix from the same group of matrices (potato, strawberry and olive)

Pesticide	100 µg/kg																													
	TOMATO										ORANGE										AVOCADO									
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G									
2,4'-DDE	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green										
2-Phenylphenol	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green										
Acephate	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red										
Acetamiprid	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	yellow										
Acrinathrin	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green										
Aldicarb	green	green	green	green	green	green	green	red	green	green	green	green	green	green	green	yellow	green	red	red	red										
Aldicarb-sulfone	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	red	red	red	red	red										

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Ametoctradin	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Azinphos-methyl	Green	Green	Green	Green	Green	Yellow	Red	Red	Red	Red	Red	Red	Red	Red	Green	Red	Green	Green	Red	Red	
Azoxystrobin	Green	Green	Green	Green	Green	Red	Red	Red	Red	Red	Red	Red	Red	Red	Green	Green	Green	Green	Green	Green	
Benalaxyl	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Benzovindiflupyr	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Yellow	Green	Green	Green	Green	Green	
Bifenthrin	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Red	Red	Green	Red	Red	Green	Green	Red	Red	
Bitertanol	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Yellow	Green	Green	Green	Green	Green	
Boscalid	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Bromopropylate	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Bupirimate	Green	Green	Green	Green	Green	Red	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	
Buprofezin	Green	Green	Green	Green	Green	Red	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Carbaryl	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	
Carbendazim	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	
Carbofuran	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	
Chlorantraniliprole	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	Green	Green	Yellow	Green	Green	Green	Yellow	
Chlorfenapyr	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Chlorfluazuron	Green	Red	Green	Green	Green	Green	Green	Yellow	Red	Yellow	Red	Red	Red	Red							
Chlорidazon	Green	Green	Green	Green	Green	Green	Green	Green	Red	Red	Red	Red	Red	Red							
Chlorpropham	Green	Green	Green	Green	Green	Red	Green	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	
Chlorpyrifos	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	
Chlorpyrifos-methyl	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Clomazone	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Clopyralid	Green	Green	Green	Green	Green	Green	Green	Yellow	Green	Green	Green	Green	Green	Green							

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Cyflufenamid																					
Cyfluthrin																					
Cyhalofop-butyl																					
Cymoxanil																					
Cypermethrin																					
Cyproconazole																					
Cyprodinil																					
Cyromazine																					
Deltamethrin																					
Diazinon																					
Dichloran																					
Dichlorvos																					
Dicofol, o, p'-																					
Dicofol, p, p'-																					
Dieldrin																					
Diethofencarb																					
Difenoconazole																					
Diflubenzuron																					
Dimethoate																					
Dimethomorph																					
Diniconazole																					
Dinotefuran																					
Diuron																					

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Emamectin B1a benzoate																					
Emamectin B1b benzoate																					
Endosulfan sulfate																					
Endosulfan-alpha																					
Endosulfan-beta																					
Epoxiconazole																					
Ethion																					
Ethirimol																					
Ethoprophos																					
Etofenprox																					
Etoxazole																					
Famoxadone																					
Fenamidone																					
Fenamiphos																					
Fenamiphos-sulfone																					
Fenamiphos-sulfoxide																					
Fenarimol																					
Fenazaquin																					
Fenbuconazole																					
Fenhexamid																					

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Fenitrothion	green	yellow	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	
Fenobucarb	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenoxy carb	green	green	green	green	green	green	green	green	green	green	green	green	yellow	green	green	green	green	green	green	green	
Fenpicoxamid	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenpropathrin	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenpropidin	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenpropimorph	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenpyrazamine	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenpyroximate	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fenthion	green	green	green	green	green	green	green	red	red	green	green	green	green	green	green	green	green	green	green	green	
Fenthion-sulfone	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	red	green	green	
Fenthion-sulfoxide	green	green	green	green	green	green	green	green	green	green	green	green	yellow	green	green	green	green	green	green	red	
Fenvalerate	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fipronil	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Flonicamid	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	green	
Fluazifop-p-butyl	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	green	red	
Flubendiamide	green	green	green	green	green	yellow	red	green	green	green	green	green	green	green	green	green	green	green	green	yellow	
Fludioxonil	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	green	green	green	
Fluensulfone	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	yellow	green	green	
Flufenacet	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	green	green	
Flufenoxuron	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	red	red	green	green	
Fluopicolide	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	
Fluopyram	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	green	

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Fluquinconazole																					
Flusilazole																					
Flutriafol																					
Fluxapyroxad																					
Formetanate																					
Hydrochloride																					
Fosthiazate																					
Haloxyfop																					
Haloxyfop-methyl																					
Heptachlor																					
Hexaconazole																					
Hexythiazox																					
Imazalil																					
Indoxacarb																					
Ioxynil																					
Iprovalicarb																					
Isocarbophos																					
Isoprothiolane																					
Isopyrazam																					
Isoxaflutole																					
Kresoxim-methyl																					
Lambda-Cyhalothrin																					

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Linuron																					
Lufenuron																					
Malathion																					
Mandipropamid																					
Mepanipyrim																					
Metaflumizone																					
Metalaxyd																					
Metconazole																					
Methamidophos																					
Methidathion																					
Methiocarb																					
Methomyl																					
Methoxyfenozide																					
Metrafenone																					
Molinate																					
Monocrotophos																					
Myclobutanil																					
Novaluron																					
Omethoate																					
Oxadiargyl																					
Oxadixyl																					
Oxamyl																					
Oxasulfuron																					

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Oxathiapipronil																					
Paclobutrazol																					
Parathion-methyl																					
Penconazole																					
Pencycuron																					
Pendimethalin																					
Penflufen																					
Penthiopyrad																					
Permethrin																					
Phenthroate																					
Phosmet																					
Picolinafen																					
Pirimicarb																					
Pirimicarb-desmethyl																					
Pirimiphos-methyl																					
Propaquizafop																					
Prothiofos																					
Pyriofenone																					
Quinalphos																					
Quinoclamine																					
Quintozene																					
Quizalofop																					

Pesticide	100 µg/kg																				
	TOMATO							ORANGE							AVOCADO						
	A	B	C	D	E	F	G	A	B	C	D	E	F	G	A	B	C	D	E	F	G
Quizalofop-ethyl	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	Green	Green	Red	Green	Green	Green	Green	
Rotenone	Green	Green	Green	Green	Green	Green	Green	Red	Red	Red	Red	Red	Red	Green	Green	Green	Green	Green	Green	Green	
Tetramethrin	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	Green	
Tolfenpyrad	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Red	Green	Green	Green	Green	Green	Green	Green	
Triallate	Yellow	Yellow	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	
Triflumizole	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	Green	
Triticonazole	Yellow	Green	Yellow	Green	Red	Red	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Yellow	Green	Green	Green	
Tritosulfuron	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Red	Green	
Zoxamide	Red	Green	Green	Green	Red	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Yellow	Yellow	Red	Green	

**A** = multiple-level calibration of the same matrix (tomato, orange and avocado); single-level calibration: **B** = 5 µg/kg, **C** = 10 µg/kg, **D** = 50 µg/kg, **E** = 100 µg/kg, **F** = 200 µg/kg; **G** = multiple-level calibration of another matrix from the same group of matrices (potato, strawberry and olive)

**Table 3.** Percentage of analysed compounds that fall into each of the classifications at 10 µg/kg (green, orange or red).

Quantification	10 µg/kg								
	TOMATO			ORANGE			AVOCADO		
	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%
Multiple-level calibration	98%	1%	1%	98%	1%	1%	93%	4%	3%
Single-level calibration: 5 µg/kg	92%	5%	3%	93%	3%	4%	89%	5%	6%
Single-level calibration: 10 µg/kg	99%	1%	0%	98%	1%	1%	86%	6%	8%
Single-level calibration: 50 µg/kg	99%	1%	0%	97%	1%	2%	93%	2%	5%
Single-level calibration: 100 µg/kg	98%	1%	1%	79%	2%	19%	92%	2%	6%
Single-level calibration: 200 µg/kg	91%	5%	4%	90%	4%	6%	91%	3%	6%
Multiple-level calibration of another matrix of the	96%	1%	3%	64%	5%	31%	64%	10%	26%

Quantification	10 µg/kg								
	TOMATO			ORANGE			AVOCADO		
	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%
same group of matrices									
Procedural calibration	97%	1%	2%	78%	8%	14%	95%	1%	4%
Sample standard addition	92%	4%	4%	86%	5%	9%	83%	4%	13%
Extract standard addition	66%	11%	23%	55%	8%	37%	45%	12%	43%

Table 4. Percentage of analysed compounds that fall into each of the classifications at 100 µg/kg (green, orange or red).

Quantification	100 µg/kg								
	TOMATO			ORANGE			AVOCADO		
	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%
Multiple-level calibration	97%	2%	1%	94%	2%	4%	98%	1%	1%
Single-level calibration: 5 µg/kg	93%	4%	3%	89%	3%	8%	88%	3%	9%

Quantification	100 µg/kg								
	TOMATO			ORANGE			AVOCADO		
	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%	0-21%	21-25%	> 25%
Single-level calibration: 10 µg/kg	98%	1%	1%	94%	1%	5%	90%	4%	6%
Single-level calibration: 50 µg/kg	97%	2%	1%	95%	0%	5%	98%	1%	1%
Single-level calibration: 100 µg/kg	99%	1%	0%	79%	6%	15%	97%	2%	1%
Single-level calibration: 200 µg/kg	95%	2%	3%	92%	4%	4%	97%	1%	3%
Multiple-level calibration of another matrix of the same group of matrices	93%	4%	3%	71%	5%	24%	64%	4%	30%

## 7. Conclusions

The evaluation of the results obtained for different matrix groups, taking into account the workload of the laboratory, leads us to the conclusion that the best way to quantify is with a **multiple-level calibration** in the same type of matrix. With this method, quantifying the **10 µg/kg** sample, **98%** of the compounds analysed in tomato and orange showed a variation with respect to the theoretical concentration of less than |21%|, and with a complicated matrix such as avocado, **93%**.

Quantification with **single-level calibration** is in general valid, obtaining in most of the cases variations above **90%**, although better quantification is obtained when the level used to quantify is **very close** in concentration to the fake sample. As this cannot be known in advance, quantification may be somewhat affected.

Quantification using a multiple-level calibration prepared with **another matrix** from the same commodity group may be acceptable in simpler matrix groups. However, in more complex matrices it seems to be less advisable, as the differences in quantification can increase. In the case of tomato quantified with the potato multiple-level calibration by GC-MS/MS, 100% of the pesticides presented a variation of less than |21%|. For the compounds analysed by LC-MS/MS, similar results were found, 94%. In the case of **avocado**, 70% of the compounds analysed by GC-MS/MS obtained values lower than |21%| quantified with olive multiple-level calibration, and for LC-MS/MS the percentage of compounds analysed below this value was 62%. The most remarkable case is **orange**, where the percentage of compounds analysed by GC-MS/MS with a percentage of variation lower than |21%| quantified with strawberry multiple-calibration was 96%, and 52% for the compounds analysed by LC-MS/MS.

Analysis by **procedural calibration** is very acceptable: 97% of the compounds analysed in tomato, 78 % in orange and 95% in avocado showed a variation of less than |21%| quantifying the 10 µg/kg level. In addition, this method of quantification can **compensate differences in recoveries**. However, quantifying with this method requires the use of significant **amounts** of analytical standards.

Quantification by **standard addition to the sample** gives generally acceptable results, obtaining variations between **92%-83%** in the case of the 10 µg/kg level, and may be a perfectly applicable option when the laboratory does not have an acceptable **blank sample** to prepare a matrix-matched calibration curve.

Results obtained in this study by **standard addition to the extract** are generally **lower**, as the percentage of compounds that showed a variation below | 21% | at the 10 µg/kg level were below **70%** in all matrices evaluated.

## 8. References

- [1] Analytical Quality Control and Method Validation Procedures for Pesticide Residues Analysis in Food and Feed. Document N° SANTE/11312/2021. EUROPEAN UNION REFERENCE LABORATORIES.
- [2] COMMISSION IMPLEMENTING REGULATION (EU) 2022/741 of 13 May 2022 concerning a coordinated multiannual control programme of the Union for 2023, 2024 and 2025 to ensure compliance with maximum residue levels of pesticides and to assess the consumer exposure to pesticide residues in and on food of plant and animal origin.
- [3] Working document on pesticides to be considered for inclusion in the national control programmes to ensure compliance with maximum residue levels of pesticides residues in and on food of plant and animal origin (SANCO/12745/2013).

## APPENDIX I: MASS TRANSITIONS

**Table A.** Detection parameters for the selected compounds analysed by LC-MS/MS.

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
Acephate	184	143	184	125	2,793
Acetamiprid	223	126	223	56	6,076
Aldicarb	213	89	213	116	7,472
Aldicarb-sulfone	239,9	223	239,9	86	3,545
Ametoctradin	276,2	176,1	276,2	149	12,978
Azinphos-methyl	318	132,1	318	261	10,331
Azoxystrobin	404	372	404	344	10,727
Benzovindiflupyr	398	377,9	398	342	12,441
Bifenthrin	440,1	181	440,1	198,2	15,151
Bitertanol	338,2	269,2	338,2	99,1	12,655
Buprofezin	306	201	306	116	13,568
Carbaryl	202	145	202	127	9,002
Carbendazim	192	160	192	132	4,213
Carbofuran	222	165	222	123	8,674
Chlorantraniliprole	483,9	452,9	483,9	285,9	10,345
Chlorfluazuron	540	382,9	540	158,1	14,224
Chloridazon	222,1	104,1	222,1	92	6
Chlorpropham	214	123,9	214	96	8,581
Clomazone	240,1	124,9	240,1	127,8	10,498
Clopyralid	192	146	192	110	3,101
Cyflufenamid	413	294,9	413	240,8	12,81
Cyhalofop-butyl	375,1	256	375,1	120,1	13,061
Cymoxanil	199,1	128	199,1	110,9	6,416
Cyprodinil	226,2	92,9	226,2	76,9	11,634
Cyromazine	167	59,9	167	125	1,761
Diazinon	305	169	305	153	12,535
Dichlorvos	220,8	108,8	220,8	78,9	8,512
Diethofencarb	268	226	268	180	10,639
Difenoconazole	406	251	406	337	12,883
Diflubenzuron	311	158	311	141	11,881
Dimethoate	230	199	230	171	6,054
Dimethomorph	388	301	388	165	11,053
Diniconazole	326,1	70	326,1	159	12,987
Dinotefuran	203,1	129,1	203,1	114,1	3,321
Diuron	233	72,1	233,03	160	10,035
Emamectin B1a benzoate	886,5	158,1	886,5	302,2	13,431
Emamectin B1b benzoate	872,5	157,9	872,5	82	13,164
Epoxiconazole	330,1	121	330,1	101,2	11,761

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
Ethion	385,1	199	385,1	171	13,713
Ethirimol	210,2	140,1	210,16	43,1	7,35
Ethoprophos	243,1	130,9	243,1	97	11,87
Etofenprox	394,2	177,3	394,2	359,1	14,887
Etoxazole	360	140,9	360	304	14,047
Famoxadone	392	331	392	238	12,486
Fenamiphos	304,1	217,1	304,1	234	12,078
Fenamiphos-sulfone	336,1	266	336,1	188	9,012
Fenamiphos-sulfoxide	320,1	108,1	320,11	292,1	8,805
Fenarimol	331	268	331	259	11,761
Fenazaquin	307,3	161,3	307,3	147,2	14,367
Fenbuconazole	337,1	70	337,1	125,1	11,917
Fenhexamid	302	97	302	55	11,645
Fenobucarb	207,8	95	207,8	152,1	10,832
Fenoxy carb	302,2	88,2	302,2	116,2	12,033
Fenpicoxamid	615,3	515	615,3	238,9	13,24
Fenpropidin	274,3	147,1	274,3	85,8	10,362
Fenpropimorph	304,3	147,1	304,3	130	10,645
Fenpyrazamine	332,2	272,1	332,2	230,2	11,508
Fenpyroximate	422,2	366,2	422,21	107	13,951
Fenthion	279	247,1	279	169,1	12,264
Fenthion-sulfone	310,7	125	310,7	108,8	9,271
Fenthion-sulfoxide	295	280	295,02	109	8,988
Fipronil	434,9	329,9	434,9	249,9	12,195
Flonicamid	230,1	202,6	230,1	173,9	4,378
Flubendiamide	680,9	254	680,9	273,9	12,396
Fludioxonil	265,9	228,9	265,9	158	11,101
Flufenacet	364,1	152	364,1	194,1	11,835
Flufenoxuron	489,1	158	489,1	140,9	13,988
Fluquinconazole	376	307,1	376	108	11,49
Flutriafol	302,1	70,1	302,1	95	9,839
Fluxapyroxad	381,9	362	381,9	342	11,244
Formetanate Hydrochloride	222,1	165,1	222,13	65,1	2,895
Fosthiazate	284	227,8	284	103,8	9,475
Haloxyfop	362,1	316,2	362,1	288,1	12,135
Haloxyfop-methyl	375,9	316	375,9	287,9	12,957
Hexaconazole	314,1	70,1	314,1	159	12,712
Hexythiazox	353,1	228,2	353,1	168,2	13,896
Imazalil	297	159	297	255	9,503
Ioxynil	369,8	126,8	369,8	214,8	10,039
Iprovalicarb	321,2	119	321,2	202,9	11,806

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
Isoprothiolane	291	230,7	291	189,1	11,22
Isoxaflutole	360	250,9	360	219,7	10,096
Linuron	249	160,1	249,02	133	10,73
Lufenuron	508,9	339	508,9	325,9	13,678
Malathion	331	127,1	331	285	11,307
Mandipropamid	412,1	328,1	412,13	356,1	11,095
Metaflumizone	505	302	505	328	13,402
Metalaxyll	280,3	220	280,3	192,4	10,073
Metconazole	320,1	70,1	320,1	125	12,663
Methamidophos	142,1	94,1	142,1	125	2,33
Methiocarb	226	169	226,1	121,1	10,903
Methomyl	163,1	106	163,1	88	4,137
Methoxyfenozide	369,3	149	369,3	133	11,518
Metrafenone	409,1	209,1	409,1	226,9	12,769
Monocrotophos	224,2	127	224,2	193,1	4,739
Myclobutanil	289,2	70,2	289,2	125,1	11,478
Novaluron	490,8	470,7	490,8	305,1	13,215
Omethoate	214,1	125	214,1	183	3,128
Oxadiargyl	341,1	222,9	341,05	150,9	12,68
Oxamyl	237	72	237	90	3,839
Oxasulfuron	407,1	150,1	407,1	209,7	8,158
Oxathiapipronil	540,2	500	540,2	522	11,228
Paclobutrazol	294,1	70,1	294,1	125,2	11,26
Penconazole	284	70	284	159	12,377
Pencycuron	329,1	125,1	329,1	89,1	12,901
Pendimethalin	282,1	212,1	282,1	194,1	13,806
Penflufen	318,1	234	318,1	141	12,327
Penthiopyrad	357,9	149	357,9	207,6	12,476
Phenthioate	321	247,1	321	79,1	12,216
Phosmet	318	160	317,99	133	10,331
Pirimicarb	239,2	72,2	239,2	182,1	7,603
Pirimicarb-desmethyl	225,1	72,1	225,1	168,1	5,294
Pirimiphos-methyl	306,2	164,2	306,2	108,2	12,583
Propaquizafop	444,1	99,9	444,1	371	13,376
Pyriofenone	366,1	183,9	366,1	209	12,787
Quinalphos	299,1	270,8	299,1	242,8	12,079
Quinoclamine	208	105,1	208	77	7,866
Quizalofop	345	299	345	254,9	11,76
Quizalofop-ethyl	373,1	271,2	373,09	255,1	13,186
Rotenone	395	213,1	395	192,1	11,875
Tolfenpyrad	384,1	197	384,1	170,9	13,47
Triallate	306	145	306,01	86	13,85

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
Triflumizole	346,1	277,8	346,1	72,9	13,243
Triticonazole	318,1	70,2	318,1	125,2	11,719
Tritosulfuron	446	195	446	145	10,515
Zoxamide	336	187	336	159	12,52

**Table B.** Detection parameters for the selected compounds analysed by GC-MS/MS.

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
2,4'-DDE	246	176	246	211	7,112
2-Phenylphenol	170	115	170	141	4,462
Acrinathrin	208	181	289	93	9,187
Benalaxyd	204	176	148	105	8,105
Boscalid	140	112	140	76	10,341
Bromopropylate	341	185	341	155	8,66
Bupirimate	273	193	273	108	7,528
Chlorfenapyr	247	227	247	200	7,667
Chlorpyrifos	314	258	314	286	6,499
Chlorpyrifos-methyl	286	271	288	93	6,077
Cyfluthrin	163	127	226	206	10,051
Cypermethrin	163	127	209	141	10,275
Cyproconazole	222	125	139	111	7,759
Deltamethrin	253	172	181	152,1	11,462
Dichloran	206	176	206	148	5,457
Dicofol, o, p'	251	139	139	111	8,37
Dicofol, p, p'	251	139	139	111	7,726
Dieldrin	279	243	345	263	7,48
Endosulfan sulfate	272	237	387	289	8,278
Endosulfan-alpha	239	204	195	160	7,335
Endosulfan-beta	207	172	240,9	205,9	7,858
Fenamidone	268	180	238	103	8,834
Fenitrothion	277	260	277	109	6,352
Fenpropathrin	265	210	265	89	8,696
Fenvalerate	167	125	125	89	10,828
Fluazifop-p-butyl	282	91	282	238	7,637
Fluensulfone	119	92	226	206	4,738
Fluopicolide	209	182	173	109	8,269
Fluopyram	173	145	223	196	6,928
Flusilazole	233	165	233	152	7,565
Heptachlor	272	237	272	143	6,155
Indoxacarb	203	134	264	148	11,48

Compound	Precursor Ion (1)	Product ion (1)	Precursor Ion (2)	Product ion (2)	Ret. Time (min)
Isocarbophos	136	108	230	212	6,633
Isopyrazam	159	139	359	303	9,478
Kresoxim-methyl	206	116	206	131	7,566
Lambda-Cyhalothrin	197	141	197	161	9,187
Mepanipyrim	222	207	222	158	7,259
Methidathion	145	85	145	58	7,132
Molinate	187	126	126	55	4,492
Oxadixyl	163	132	163	117	7,984
Parathion-methyl	263	109	233	124	6,124
Permethrin	163	127	183	153	9,666
Picolinafen	376	238	238	145	8,709
Prothiofos	309	239	309	221	7,33
Quintozene	295	237	295	265	5,589
Tetramethrin	164	107	164	77	8,645