

Analysis of Quaternary Ammonium Compounds (QACs) in Fruits and Vegetables using QuEChERS and LC-MS/MS

Version 5 (last update: 24.03.2016)

Brief description:

A method is presented for the analysis of Quaternary ammonium compounds (QACs) in food of plant origin. The compounds are extracted using the QuEChERS method and analysed by LC-MS/MS in the ESI (pos.) mode using a polar endcapped column and a slightly acidic gradient.

Compound details:

Quaternary ammonium compounds (QACs) are surface-active substances containing a quaternary cationic nitrogen atom, substituted by alkyl chains of varying length. They are used as biocides, pesticides, disinfectants and additives for technical applications.


BAC is a mixture of alkylbenzyltrimethylammonium chlorides of various even-numbered alkyl chain lengths (C8-C18). The greatest biocide activity is associated with the C12-C14 derivatives, which are the main components of the mixture.

DDAC is a mixture of alkyl-quaternary ammonium salts with typical alkyl chain lengths of C8, C10 and C12. C10 is the main component and makes more than 90 % of the mixture. In most cases the term DDAC is thus used for the Didecyltrimethylammonium chloride congener.

DDAC used to be approved in the EU as a plant protection product (for ornamental plants). BAC has never been approved as active substance in plant protection products. Still, as both Both BAC and DDAC have pesticide properties they were until a few years ago contained in various products used for the treatment of various food-related crops in the field within the EU. Some of these products adulterated were even marketed as suitable for organic agriculture (e.g. grapefruit seed extract based products containing >20% BAC) without the QACs contained being declared.


BAC and DDAC are furthermore widely used as **biocides** for the **cleaning/sanitation of surfaces** in the areas of food, feed and medical production sectors as well as domestically. Personal care products such as hand sanitizers, wet wipes and eye-drops may also contain them.

Dairy Industry is probably the most important food-related field of DDAC and BAC use. DDAC and BAC containing products are for example employed for the disinfection of all kind of surfaces milking equipment, milk storage tanks and ice-cream machines. BAC and DDAC are popular among dairy farmers, where they are used to disinfect udder and milking equipment with the goal of preventing mastitis and producing raw milk with low bacterial count, which is preferential both from the hygienic and economic point of view. Unlike many chlorine-containing products BAC and DDAC do not cause skin irritation.




Compound facts at a glance:

| Benzalkonium chloride (BAC) CAS 8001-54-5 | |
|---|---|
| Parameter | Value |
| Molecular Formula: | C ₉ H ₁₃ CINR (R=C ₈ H ₁₇ to C ₁₈ H ₃₇) |
| Molecular weight: | variable |
| Residue definition EU | Benzalkonium chloride (mixture of alkylbenzyl-dimethylammonium chlorides with alkyl chain lengths of C ₈ , C ₁₀ , C ₁₂ , C ₁₄ , C ₁₆ and C ₁₈) |



n=8, 10, 12, 14, 16, 18

| Didecyl-dimethylammonium chloride (DDAC-C10) CAS 7173-51-5 | |
|--|--|
| Parameter | Value |
| Molecular Formula: | C ₂₂ H ₄₈ CIN |
| Molecular weight: | 362,08 |
| Residue definition EU: | Didecyl-dimethylammonium chloride (mixture of alkyl-quaternary ammonium salts with alkyl chain lengths of C ₈ , C ₁₀ and C ₁₂) |




Legal Aspects and Monitoring Programs: In October 2012 the Standing Committee on the Food Chain and Animal Health (SCoFAH) endorsed Guidelines on measures to be taken as regards the presence of DDAC¹ and BAC² in or on food and feed. These guidelines contained an agreed temporary enforcement level of 0.5 mg/kg for food and feed. It was further recommended that EU Member States carry out investigations on the possible causes of BAC / DDAC contamination and a monitoring program was put in place to get an overview of the BAC and DDAC levels in food and feed of plant and animal origin. Based on the monitoring data collected (summarized in an EFSA report³) a flat MRL level of 0.1 mg/kg was established in 2014 for both didecyl-dimethylammonium chloride and benzalkonium chloride in all food products (Reg. 1119/2014/EU). These MRLs apply to the chloride salts. A revision of these MRLs is planned to take place by 31 December 2019 and if necessary they will be amended modified based on newer residue data. Before the agreement on temporary MRLs (0.5 mg/kg) in 2012 and the implementation of specific MRLs in 2014 (0.1 mg/kg), the default MRL of 0.01 mg/kg laid down in Reg. 396/2005/EC used to apply.

¹ European Commission, 2012a. Guidelines as regards measures to be taken as regards the presence of Didecyl Dimethyl Ammonium Chloride (DDAC) in or on food and feed agreed by the Standing Committee of the Food Chain and Animal Health (SCoFAH) on 13 July 2012 and modified on 5 October 2012

² European Commission, 2012b. Guidelines as regards measures to be taken as regards the presence of Benzalkonium Chloride (BAC) in or on food and feed agreed by the Standing Committee of the Food Chain and Animal Health (SCoFAH) on 25 July 2012 and modified on 5 October 2012

³ European Food Safety Authority; Evaluation of monitoring data on residues of didecyl-dimethylammonium chloride (DDAC) and benzalkonium chloride (BAC). EFSA supporting publication 2013:EN-483, 30 pp



Standard Materials Supply:

Table 1: Exemplary Providers of QACs (11.03.2016)⁴


| Acronym | Compound | CAS | Company | Order No. |
|-------------------------------------|---|-----------|---------|---------------|
| Native standards | | | | |
| BAC-C8 | Benzyltrimethyloctylammonium chloride | 959-55-7 | HPC | 675698 |
| BAC-C10 | Benzyltrimethyldecylammonium chloride | 965-32-2 | Aldrich | 13371 |
| BAC-C12 | Benzyltrimethyldodecylammonium chloride | 139-07-1 | Aldrich | 13380 |
| BAC-C14 | Benzyltrimethyltetradecylammonium chloride | 139-08-2 | Fluka | 13401 |
| BAC-C16 | Benzyltrimethylhexadecylammonium chloride | 122-18-9 | Sigma | B4136 |
| BAC-C18 | Benzyltrimethyloctadecylammonium chloride | 122-19-0 | HPC | 674644 |
| DDAC-C8 | Dimethyldioctylammonium bromide | 3026-69-5 | HPC | 675696 |
| | Dimethyldioctylammonium chloride | 5538-94-3 | LGC | DRE-C12726470 |
| DDAC-C10 | Didecyldimethylammonium chloride | 7173-51-5 | LGC | DRE-C12588000 |
| DDAC-C12 | Didodecyldimethylammonium chloride | 3401-74-9 | Chemos | 135628 |
| | | | TRC | D441760 |
| Internal Standards (Options) | | | | |
| BAC-C10 D7 | D7-Benzyltrimethyldecylammonium chloride | | HPC | 674610 |
| BAC-C12 D6 | D6-Benzyltrimethyldodecylammonium iodide | | HPC | 674572 |
| BAC-C14 D7 | D7-Benzyltrimethyltetradecylammonium chloride | | HPC | 674611 |
| DDAC-C10 D6 | D6-Didecyldimethylammonium iodide | | HPC | 674541 |

Measurement


Table 2: LC details for QACs (exemplary)

| | | | |
|---------------------|---|---------------|-------------|
| LC | Agilent 1290 Infinity UPLC | | |
| MS/MS | ABSCIEX API 4000 Q-Trap, run in ESI positive mode | | |
| Column | Phenomenex Aqua 5µ C18 125A | | |
| Pre-column | Aqua C18 125A 4mm x 2mm (Phenomenex AJO-7510) | | |
| Mobile Phase | A: 5 mmol NH ₄ -formiate in H ₂ O + 0,01 % formic acid B: 5 mmol NH ₄ -formiate in methanol (MeOH) + 0,01 % formic acid | | |
| Ionisation mode | ESI pos | | |
| Gradient | %A | Flow [mL/min] | Time [min] |
| | 95 | 0.4 | 0 |
| | 50 | 0.4 | 2 |
| | 40 | 0.4 | 8 |
| | 10 | 0.4 | 12 |
| | 10 | 0.4 | 14 |
| | 95 | 0.4 | 14.5 |
| | 95 | 0.4 | 21 |
| Injection volume | 5 µL | | |
| Dilution of extract | no | | |
| Column temperature | 40°C | | |

⁴ **Disclaimer:** this information is given solely for the convenience and information of the readers and it does not constitute any advertisement, preference or endorsement of the companies and products listed


Table 3: MS/MS details for QACs (ESI-negative mode, Tune-data ABSciex 4000Q) (exemplary)

| Compound | Sensitivity Ranking (1= best) | Parent Mass | Daughter Mass | DP | CE | CXP |
|-------------------------------------|-------------------------------|-------------|---------------|-----|-----|-----|
| BAC-C8 | 1 | 248 | 156 | 36 | 25 | 25 |
| | 2 | 248 | 91 | 36 | 41 | 8 |
| BAC-C10 | 1 | 276 | 184 | 55 | 27 | 10 |
| | 2 | 276 | 91 | 55 | 37 | 36 |
| BAC-C12 | 1 | 304 | 212 | 91 | 29 | 10 |
| | 2 | 304 | 91 | 91 | 37 | 16 |
| BAC-C14 | 1 | 332 | 240 | 83 | 31 | 10 |
| | 2 | 332 | 91 | 83 | 59 | 8 |
| BAC-C16 | 1 | 360 | 268 | 78 | 33 | 12 |
| | 2 | 360 | 91 | 78 | 67 | 10 |
| BAC-C18 | 1 | 388 | 296 | 31 | 33 | 8 |
| | 2 | 388 | 91 | 31 | 111 | 26 |
| DDAC-C8 | 1 | 270 | 158 | 46 | 35 | 6 |
| | 2 | 270 | 43 | 46 | 61 | 2 |
| DDAC-C10 | 1 | 326 | 186 | 61 | 39 | 12 |
| | 2 | 326 | 41 | 61 | 93 | 6 |
| DDAC-C12 | 1 | 382 | 214 | 121 | 43 | 6 |
| | 2 | 382 | 58 | 121 | 67 | 4 |
| Internal Standards (Options) | | | | | | |
| BAC-C10 D7 | - | 283 | 98 | 81 | 51 | 4 |
| BAC-C12 D6 | - | 310 | 218 | 91 | 37 | 16 |
| BAC-C14 D7 | - | 339 | 98 | 96 | 61 | 4 |
| DDAC-C10 D6 | - | 332 | 192 | 96 | 41 | 10 |
| Chlorpyrifos D10 | | 360 | 199 | 36 | 27 | 4 |



Validation Data

Table 4: Compilation of exemplary recovery data for various commodity types and spiking levels

| Commodity group | Commodity | Compound | spiking Level | Average Recovery | RSD _R (Reproducibility) | No. of recoveries |
|----------------------|-----------|----------|---------------|------------------|---------------------------------------|-------------------|
| High water content | Asparagus | BAC-C8 | 0.025 | 102 | 13 | 7 |
| | | BAC-C10 | 0.025 | 102 | 12 | 10 |
| | | BAC-C12 | 0.025 | 101 | 13 | 10 |
| | | BAC-C14 | 0.025 | 99 | 11 | 9 |
| | | BAC-C16 | 0.025 | 93 | 14 | 8 |
| | | BAC-C18 | 0.025 | 90 | 18 | 8 |
| | | DDAC-C10 | 0.025 | 98 | 13 | 10 |
| | Cucumbers | BAC-C8 | 0.025 | 100 | 8 | 11 |
| | | BAC-C10 | 0.025 | 100 | 7 | 11 |
| | | BAC-C12 | 0.025 | 100 | 10 | 11 |
| | | BAC-C14 | 0.025 | 97 | 9 | 11 |
| | | BAC-C16 | 0.025 | 95 | 8 | 11 |
| | | BAC-C18 | 0.025 | 91 | 9 | 11 |
| | | DDAC-C10 | 0.025 | 97 | 8 | 11 |
| High acid content | Mandarins | BAC-C8 | 0.025 | 100 | 4 | 5 |
| | | BAC-C10 | 0.025 | 101 | 3 | 5 |
| | | BAC-C12 | 0.025 | 103 | 5 | 5 |
| | | BAC-C14 | 0.025 | 101 | 2 | 5 |
| | | BAC-C16 | 0.025 | 102 | 3 | 5 |
| | | BAC-C18 | 0.025 | 97 | 5 | 5 |
| | | DDAC-C10 | 0.025 | 99 | 2 | 5 |
| | Oranges | BAC-C8 | 0.025 | 100 | 8 | 8 |
| | | BAC-C10 | 0.025 | 102 | 8 | 8 |
| | | BAC-C12 | 0.025 | 104 | 10 | 8 |
| | | BAC-C14 | 0.025 | 104 | 12 | 8 |
| | | BAC-C16 | 0.025 | 100 | 9 | 8 |
| | | BAC-C18 | 0.025 | 97 | 11 | 8 |
| | | DDAC-C10 | 0.025 | 103 | 4 | 6 |
| Dry (pulses cereals) | Wheat | BAC-C8 | 0.05 | 107 | 7 | 7 |
| | | BAC-C10 | 0.05 | 104 | 7 | 7 |
| | | BAC-C12 | 0.05 | 99 | 10 | 7 |
| | | BAC-C14 | 0.05 | 96 | 9 | 7 |
| | | BAC-C16 | 0.05 | 86 | 10 | 6 |
| | | BAC-C18 | 0.05 | 80 | 16 | 7 |
| | | DDAC-C10 | 0.05 | 100 | 10 | 7 |

| Commodity group | Commodity | Compound | spiking Level | Average Recovery | RSD _R (Reproducibility) | No. of recoveries |
|--------------------|-----------|----------|---------------|------------------|---------------------------------------|-------------------|
| High sugar Content | Honey | BAC-C8 | 0.083 | 103 | 5 | 6 |
| | | BAC-C10 | 0.083 | 104 | 4 | 6 |
| | | BAC-C12 | 0.083 | 105 | 5 | 6 |
| | | BAC-C14 | 0.083 | 104 | 4 | 6 |
| | | BAC-C16 | 0.083 | 103 | 4 | 6 |
| | | BAC-C18 | 0.083 | 103 | 4 | 6 |
| | | DDAC-C10 | 0.083 | 103 | 4 | 6 |

The sensitivity of detection is good so that determinations below 0.01 mg/kg are principally possible. It should be considered, however, that the compounds are ubiquitous and have the tendency to strongly interact with surfaces and thus the potential for cross-contaminations. Background levels mostly concern BAC-C12 and BAC-C14, which are typically the predominant components in formulations. Thus LOQs of 0.02 mg/kg are more realistic.

Recoveries can also be found at www.eurl-pesticides-datapool.eu

Document History

| Date | Action | Changes |
|-------------------|-------------------|---|
| 29.06.2012 | Publication of V1 | |
| 05.07.2012 | Publication of V2 | Addition of detailed LC-Conditions |
| 06.07.2012 | Publication of V3 | Elimination of errors |
| 14.09.2012 | Publication of V4 | |
| 24.03.2016 | Publication of V5 | General text revision Revision of legal background information Inclusion of BAC-C8, BAC-C18, DDAC-C8, DDAC-C12, and isotope labelled internal standards Update of LC and MS/MS method details Update of recovery data |