

## EUPT-FV-SM13 SPECIFIC PROTOCOL

### European Union Proficiency Test for Pesticide Residues in Fruits and Vegetables Screening Multiresidue Methods (2021)

#### Introduction

This protocol is complementary to the General Protocol for EU Proficiency Tests (EUPT) dealing with Pesticide Residues in Food and Feed. This Proficiency Test is organised by the EURL for Pesticide Residues in Fruits and Vegetables and covers the screening of pesticides using multiresidue methods of analysis.

The aim of this test is to evaluate laboratory capability when using large-scope quantitative and/or screening methods during routine analysis, for detecting and identifying unexpected pesticides at levels at, or above 0.01 mg/kg – included in and/or in addition to the laboratories' quantitative methods used for frequently detected pesticides. A second aim is to encourage official laboratories to extend the scope of their methods in a cost-effective way, by using the different MS instruments/software and methods available (whether old or new).

The evaluation of this PT will be based on qualitative information, although an **estimated concentration** will be requested for those pesticides that are detected, **only for informative purposes**.

#### Test item

This proficiency test is based on the pesticide-residue analysis of **AUBERGINE**.

The pesticide treatments will be carried out post-harvest using either commercial formulation in micro-spray solutions or using standard solutions. The test item will be frozen (using liquid nitrogen), chopped, homogenised and sub-sampled into polyethylene bottles that have previously been coded.

Ten of these bottles containing the test item will be chosen randomly and analysed to check for homogeneity.

The test item will be stored frozen (-20°C) prior to shipment to participants.

The stability tests will be carried out by the EURL-FV laboratory at the University of Almería (accredited under ISO/IEC 17025 by the Spanish accreditation body, ENAC).

Shortly before the test item shipment, three bottles that are stored in the freezer at -20°C will be chosen randomly and stored in a -80°C freezer. After the deadline for reporting

results, those three bottles stored at  $-80^{\circ}\text{C}$ , together with three other bottles that are stored in the freezer at  $-20^{\circ}\text{C}$  and will be chosen randomly will be analysed by duplicate.

Moreover, regarding the stability of the sample arriving not completely frozen, a duplicate analysis of three bottles reproducing the delivery conditions that the samples experienced for 48 hours will be performed. Laboratories could therefore be sufficiently confident in accepting the treated test item even if it was not completely frozen.

These results will not be included in the statistical analysis of the proficiency test. The aim is solely to check pesticide stability during the shipping process and for the duration of the proficiency test.

### Steps to follow

This Proficiency Test will be made up of the following six essential steps:

1. To participate, each laboratory must complete the Application Form online, available on the EURL-FV Web page, before the deadline stipulated on the calendar.
2. Laboratories will then receive an e-mail confirming their participation in this exercise, assigning them a unique Laboratory Code. With this code, laboratories will fill in an excel file including their results.
3. The sample delivery will be free of charge to those laboratories already participating in **EUPT-FV23**. For those who are not **EUPT-FV23** participants, please see **cost for shipment of the test item** for further details. **Payments without an Invoice Number to identify them will not be considered paid.**
4. When the participant laboratories receive the test item (and not before), they must fill in the "sample reception form" or send an e-mail to [omalato@ual.es](mailto:omalato@ual.es) to inform the Organiser that they have accepted the test item. If no test item has been received by **9<sup>th</sup> February 2021**, please contact the Organiser by e-mail ([cferrer@ual.es](mailto:cferrer@ual.es) and [omalato@ual.es](mailto:omalato@ual.es)).
5. The participating laboratories must respect the deadline for submitting the results. Results must be reported using an excel file: **Form 1 - Results** within 72 hours after the arrival of the test item, and this excel file must be sent to [omalato@ual.es](mailto:omalato@ual.es) before the 72 h deadline. If this deadline were during the weekend, it would be postponed to the next working day.
6. The Organiser will evaluate the results at the end of the proficiency test, once the deadline for the receipt of results has passed. The Organiser will prepare a Preliminary Report that will be sent to the participants and uploaded to the website

to show the pesticides reported. After the revision of all the data by the Scientific Panel, a Final Report will be prepared, and the organiser will upload an electronic version on the EURL-FV website. This report will include information regarding the design of the test, the homogeneity and stability test results, an evaluation of the participant's results as well as graphical displays of the results and any conclusions. Any other relevant information considered of value may also be included.

### Amount of Test Item

Participants will receive:

- Approximately 200 g of **AUBERGINE** test item treated with pesticides.

### Shipment of Test Items

All test items will be frozen and packed in polystyrene boxes surrounded by dry ice and packed into cardboard boxes.

The shipment of the test items will start on **1<sup>st</sup> February 2021**. Laboratories must make their own arrangements for the receipt of the package. They must inform the Organiser of any public holidays in their country/city during the delivery period given in the calendar as well as making the necessary arrangements to receive the shipment, even if the laboratory is closed.

### Advice on Test Item Handling

Once received, the test item should be stored deeply frozen (-18°C or less) prior to analysis to avoid any possible deterioration/spoilage. The test item should be mixed thoroughly before taking the analytical portion(s).

All participants should use their own routine standard operating procedures for extraction, clean-up and analytical measurement and their own reference standards for identification.

### Test Item Receipt

Once the laboratory has received the test item, the Organiser must be notified filling in the "sample reception form" available in our test webpage, or by e-mail to [omalato@ual.es](mailto:omalato@ual.es) indicating the date of receipt, and acknowledging its acceptance. If the laboratory does not inform the Organiser by 10<sup>th</sup> February 2021 (at the latest) stating that no sample has been received, the Organiser will assume that the test item has been received and accepted.

**Form 1 – Results**

The evaluation of this PT will be based on qualitative information, although an estimated concentration will be requested for those pesticides that are detected, only for informative purposes.

It has been decided by the Quality Control Group, and based on the received questionnaires, that a target pesticide list will not be provided.

Laboratories must fill in their results in “Results Form” (an excel file).

On this form, the laboratory should report the name of each of the pesticides detected. Each pesticide may be reported more than once if it has been detected by more than one method or identification criteria, as long as details of each method used are also provided.

Information on the parameters and/or criteria used for detecting and reporting the pesticides found will be requested, such as deviation from expected retention time, and MS identification details.

The idea is to ascertain if the methods are used in routine or just specifically for this test and if the identification is undertaken manually or automatically. Moreover, the range over which the method operates will be required (the concentration range, - the minimum and maximum level of the screening method that is used to detect pesticides)

Information **must be sent up to 72 hours after sample arrival** in the laboratory. After the deadline, results submission and/or changes to the results form will no longer be possible.

**EUPT-FV-SM13 CALENDAR**

Activity	Date
Receiving Application Form from invited laboratories.	4 <sup>th</sup> Dec 2020 – 8 <sup>th</sup> Jan 2021
Specific Protocol published on the Web site.	20 <sup>th</sup> Jan 2021 at the latest
Sample distribution.	1 <sup>st</sup> February 2021
Deadline for receiving results: Fill in “Results form	72 hours after receiving the sample
Preliminary Report	Last week of March 2021
Final Report distributed to the Laboratories.	December 2021*
Receiving Application Form from invited laboratories.	4 <sup>th</sup> Dec 2020 – 8 <sup>th</sup> Jan 2021

### Cost for shipment of the test item

Only those laboratories not participating in **EUPT-FV22** will have to pay the fee for sample shipment: 250 euros for EU/EFTA laboratories and 350 euros for other participants. For the payment procedures, each laboratory can specify their details and requests for invoices when applying for the test.

**Please, do not pay for this EUPT until we send you the invoice.**

**Remember to include your Invoice number in the subject of the bank transfer.**

Payment details are as follows:

**BANK NAME: CAJAMAR - Caja Rural Sociedad Corporativa de Crédito**  
**BANK ACCOUNT HOLDER: Universidad de Almería**  
**BANK ADDRESS: Office Number 990. Universidad de Almería. Spain**  
**ACCOUNT NUMBER: ES0730580130172731005000**  
**SWIFT: CCRIES2A**  
**REFERENCE GIVEN: INVOICE NUMBER**

### Contact information

The official organising group details are as follows:

Universidad de Almería. Edificio Químicas CITE I  
Ctra. Sacramento s/n  
04120, La Cañada de San Urbano  
Almería - Spain  
Fax No.: +34 950015008

### Organising team (e-mail and phone no.):

Dr. Amadeo R. Fernández-Alba	EURL-FV amadeo@ual.es	+34 950015034
Dr. Carmen Ferrer Amate	EURL-FV cferrer@ual.es	+34 950014102
Mr. Octavio Malato Rodríguez	EURL-FV omalato@ual.es	+34 950214423
Mr. Francisco José Díaz Galiano	EURL-FV diaz-galiano@ual.es	+34 950015645

### Quality Control Group

Dr. Antonio Valverde, University of Almería, Spain  
Dr. Paula Medina, European Food Safety Authority, Parma, Italy.

### Advisory Group

Dr. Michelangelo Anastassiades, CVUA, Stuttgart, Germany.  
Dr. Magnus Jezussek, Bavarian Health and Food Safety Authority, Erlangen, Germany.  
Dr. André de Kok, Wageningen Food Safety Research, Wageningen, The Netherlands.  
Mr. Ralf Lippold, CVUA, Freiburg, Germany.  
Dr. Sonja Masselter, AGES GmbH Institute for Food Safety, Innsbruck, Austria.  
Dr. Hans Mol, Wageningen Food Safety Research, Wageningen, The Netherlands.  
Mr. Finbarr O'Regan, Pesticide Control Laboratory, Celbridge, Ireland.  
Ms. Patrizia Pelosi, Istituto Superiore di Sanità, Roma, Italy.  
Dr. Tuija Pihlström, National Food Agency, Uppsala, Sweden.  
Dr. Mette Erecius Poulsen, DTU, Copenhagen, Denmark.