

# Ongoing method Performance verification – Procedural calibration

Hermann Unterluggauer

**AGES - Institute for Food Safety, Innsbruck**



O-14 @"Themed day"

# FATchers – a novel, simple and fast sample preparation in food of animal origin

Hermann Unterluggauer

**AGES - Institute for Food Safety, Innsbruck**

# In-house-method „FATchers“

Procedural Standard Calibration

## Extraction

- **5g sample \***
- **add 20ml mix ACN:H2O (1:1) + ISTD – shake**
- **add 5ml n-Hexan – shake**
- **add Quechers salt 1 #**  
– shake, centrifuge

# 450 mg MgSO<sub>4</sub>, 75 mg CHROMABOND® Diamino



## Clean up

- **dispersive SPE - add Quechers salt 2**  
(600 mg MgSO<sub>4</sub>, 300 mg CHROMABOND® Diamino, 450 mg CHROMABOND® C18ec)
- – 2 min shake (Collomix) and centrifuge



## Analysis

- GC-MS/MS (GC-MSD NCI) and/or LC-MS/MS

*Optional:*

- **Take 2 mL of the eluate and evaporate** to dryness
- **Reconstitute in 250µl ACN (LC) or 250µl Cyclohexane (GC) and inject**

\* modified method for honey: 10g of sample weight - resolve honey in 10ml water, add 10ml ACN)

# FATchers – more than a compromise



Suitable for quite an extended range of **F**at **A**nd **T**ricky matrices



© Aleksey Mnogosmyslov\_Shutterstock.com



© LumenSt\_Shutterstock.com

... plant oils



© Alexander Prokopenko/Shutterstock.com

insects



© Lightspring\_Shutterstock.com

# Administration of QC-Samples in AGES Innsbruck

**Procedural calibration** is an integral part of FATchers method as it **compensates for losses due to liquid/liquid partitioning step**

- Matrix-matched calibration levels prepared by spiking of sample aliquots before sample extraction
- ACN-Extracts are stored in ALS vials over several month
- PC calibration levels are used for each batch of samples, together with at least one freshly prepared QC sample at LOQ-level



# Overview – Data set I



<b>GC-MS/MS</b>	<b>plant oils</b>	<b>honey</b>	<b>milk</b>
spike level [mg/kg]	0,02	0,01	0,01
QC sample number(on average)	30	105	45
total number of compounds	370	380	370
compounds within mean apparent recovery 80 - 120 %, RSD ≤20%	99% (79 - 107%)	100% (90 - 111)	99% (77 - 108%)

E4 As a practical approach, residues results do not have to be adjusted for method bias when the mean bias is less than 20% and the default expanded measurement uncertainty of 50% is not exceeded.

# Details – Data set I

Percent of QC-samples with apparent recovery

80-120% / < 80% / > 120%



GC-MS/MS	plant oils (30 samples)	honey (105 samples)	milk (45 samples)
<b>Organochlorine</b>			
cis-chlordane	87 / 10 / 3	94 / 1 / 5	92 / 2 / 6
trans-chlordane	94 / 6 / 0	89 / 2 / 9	92 / 3 / 5
endosulfansulfat	88 / 9 / 3	92 / 8 / 0	80 / 6 / 14
gamma-HCH (Lindane)	93 / 7 / 0	82 / 7 / 11	86 / 0 / 14
Hexachlorobenzene (HCB)	100 / 0 / 0	96 / 2 / 2	98 / 2 / 0
<b>Organophosphorous compounds</b>			
Chlorpyrifos	94 / 3 / 3	95 / 3 / 2	97 / 0 / 3
Diazinon	91 / 3 / 6	86 / 3 / 11	84 / 3 / 13
Methidathion	87 / 13 / 0	84 / 12 / 4	88 / 6 / 6
Pyrazophos	90 / 10 / 0	81 / 17 / 2	83 / 8 / 9
<b>Pyrethroids</b>			
Bifenthrin	79 / 21 / 0	89 / 9 / 2	85 / 6 / 9
Cypermethrin	79 / 15 / 6	86 / 10 / 4	88 / 3 / 9
Deltamethrin	91 / 9 / 0	80 / 16 / 4	78 / 14 / 8
Permethrin	88 / 8 / 4	80 / 15 / 5	79 / 13 / 8
lambda-Cyhalothrin	84 / 16 / 0	86 / 11 / 3	83 / 10 / 7
<b>Polar (LC-amenable) compounds</b>			
Boscalid	84 / 13 / 3	88 / 10 / 2	86 / 6 / 8
Malathion	88 / 6 / 6	93 / 6 / 1	87 / 6 / 7
Phosmet	80 / 16 / 4	75 / 23 / 2	84 / 14 / 2
Trifloxystrobin	83 / 7 / 10	79 / 14 / 7	92 / 0 / 8

# Overview – Data set II



<b>GC-MS/MS</b>	<b>muscle</b>	<b>fish</b>	<b>hens egg</b>
spike level [mg/kg]	0,01	0,01	0,02
QC sample number(on average)	60	45	25
total number of compounds	370	370	370
compounds within mean apparent recovery 80 - 120 % RSD ≤20%	100% (91 - 112)	100% (87 - 115)	99% (73 - 111%)



# Details – Data set II

Percent of QC-samples with apparent recovery

80-120% / < 80% / >120%



GC-MS/MS	muscle (60 samples)	fish (45 samples)	hens egg (25 samples)
<b>Organochlorine</b>			
cis-chlordane	89 / 9 / 2	89 / 2 / 9	93 / 0 / 7
trans-chlordane	89 / 9 / 2	77 / 15 / 8	85 / 0 / 15
endosulfansulfat	65 / 19 / 16	77 / 16 / 7	48 / 38 / 14
gamma-HCH (Lindane)	79 / 12 / 9	83 / 7 / 10	85 / 4 / 11
Hexachlorobenzene (HCB)	98 / 0 / 2	90 / 8 / 2	100 / 0 / 0
<b>Organophosphorous compounds</b>			
Chlorpyrifos	94 / 1 / 5	85 / 15 / 0	85 / 0 / 15
Diazinon	86 / 5 / 9	86 / 10 / 4	88 / 0 / 12
Methidathion	79 / 8 / 13	68 / 21 / 11	100 / 0 / 0
Pyrazophos	81 / 12 / 7	76 / 14 / 10	92 / 4 / 4
<b>Pyrethroids</b>			
Bifenthrin	92 / 3 / 5	77 / 6 / 17	96 / 0 / 4
Cypermethrin	82 / 7 / 11	75 / 22 / 3	81 / 8 / 11
Deltamethrin	75 / 16 / 9	71 / 25 / 4	83 / 4 / 13
Permethrin	80 / 13 / 7	70 / 16 / 14	96 / 0 / 4
lambda-Cyhalothrin	67 / 15 / 18	80 / 13 / 7	78 / 9 / 13
<b>Polar (LC-amenable) compounds</b>			
Boscalid	92 / 5 / 3	78 / 4 / 18	92 / 4 / 4
Malathion	80 / 10 / 10	77 / 16 / 7	100 / 0 / 0
Phosmet	69 / 22 / 9	63 / 11 / 26	83 / 13 / 4
Trifloxystrobin	93 / 2 / 5	86 / 0 / 14	88 / 0 / 12

Österreichische Agentur für Gesundheit  
und Ernährungssicherheit GmbH



**Hermann Unterluggauer**

---

Technikerstrasse 70  
A - 6020 Innsbruck

**[www.ages.at](http://www.ages.at)**

Copyright © 2023 AGES/ Dr. Hermann Unterluggauer

Alle Rechte vorbehalten. Die Inhalte sind geistiges Eigentum der AGES. Diese dürfen ausschließlich für den privaten Gebrauch verwendet werden. Alle anderen Werknutzungsarten, einschließlich der Vornahme von Änderungen und Bearbeitungen, sowie eine Weitergabe an Dritte sind untersagt.